Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Made in Romania

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Technical Information Silicon based precision condensation curing impression material.

Instructions for use (GB)

Optosil®/ Xantopren® Comfort®	EN ISO 4823	Colours	Dosage		Mixing-time	Working-time ¹	Time in mouth
			Activator universal Plus				
Optosil® Comfort® Putty	Type 0	yellow	1 spoon (10,7 g/9,5 ml)	6-8 drops (0,14 g/0,13 ml)	6-8 gradations ² (0,32 g/0,29 ml)	30	1:00
Xantopren® Comfort® Medium	Type 3	green	/	/	/	/	2:30
Xantopren® Comfort® Light	Type 3	blue	/	/	/	/	3:00

*1 At room temperature of 23°C (73°F), 50% rel. humidity.

Higher temperatures shorten, lower temperatures prolong these times.

*2 Graduation marks as indicated on the mixing pad and mixing spatula (1 graduation = 1,3 cm)

The stated dosage should be observed.

Contraindications: The use of this product is contraindicated in case of known or suspected allergies against components of this product. Optosil Comfort contains peppermint oil.

Side effects: This product or one of its components may in particular cases cause hypersensitivity reactions. If suspected, information on the ingredients can be obtained from the manufacturer. Surplus material on the impression tray or in the patient's mouth should be avoided as this may result in gagging, swallowing or aspiration. To further reduce these risks during impression taking, the patient should be positioned as upright as possible.

User Advice

The Optosil®/Xantopren® assortment is a complete system.

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