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ZAKK[®] Implant

CE 0482

Hersteller

R-dental Dentalerzeugnisse GmbH - Winterhuder Weg 88 - D-22085 Hamburg
 Telefon +49 (0)40 - 22 75 76 17 - Fax +49 (0)40 - 22 75 76 18 - www.r-dental.com

ZAKK[®] Implant (KwikMix[®]-System 4:1)
Bestell-Nr.: ZKS2080

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ZAKK[®] Implant (KwikMix[®]-System 4:1)
Order-No.: ZKS2080

Description and indication

ZAKK[®] Implant in a 4:1-double-chamber-syringe (KwikMix[®]) is an universally applicable, automixing, non-eugenol, self-adhesive, semi-permanent luting cement on methacrylate basis for cementing crown and bridges on implants.

Characteristics

ZAKK[®] Implant can be light-cured and is self-curing at the same time (dual-cure).

ZAKK[®] Implant expands slightly, which results in a good seal between abutment and supraconstruction. After polymerization ZAKK[®] Implant shows an elastic behaviour, which is a great advantage in case of an optional revision.

ZAKK[®] Implant guarantees an excellent adhesion and removability.

ZAKK[®] Implant polymerizes to an elastic semi-permanent cement. The low layer thickness between abutment and supraconstruction enables a fitting position without micromovements.

Recommended procedures**1. Preparing the double-chamber-syringe**

First remove the original cap (do not re-use!) of the double-chamber-syringe and as a check press out a small amount of material without mixing cannula. Then screw on the MIXING CANNULA MPN 2 (Order-No.: MIK1046). The double-chamber-syringe is ready for use.

2. Application

Dispose first a small amount of ZAKK[®] Implant after extruding. Then introduce the material with help of the mixing cannula onto the abutment or within the supraconstruction and set it up with a slight pressure within 2 min (working time).

The intra-oral hardening time is more than 4 min (hardening time). Because of its dual-cure character ZAKK[®] Implant can be light-cured for about 20 s at light sensitive positions.

The polymerized cement can be easily separated with a scaler and alcohol.

Times mentioned for 73 °F. Higher temperatures shortens, lower temperatures prolongs these times.

Special storage and/or handling conditions

The filled MIXING CANNULA functions as a seal until the next application. Do not re-use the double chamber syringe cap. Closing the double chamber syringe by turning on the original cap between use could cause cross contamination of basis and catalyst and thus premature setting.

Important: Check the openings of the double chamber syringe and clean if necessary before installing a new mixing cannula.

Expiration date see label. Do not use after expiration date.

Do not store above 20 °C (68 °F)! Avoid exposure to sunlight.

The double chamber syringe should be stored horizontally (laying).

Precautions

Avoid contact of ZAKK[®] Implant with soft tissues, skin or eyes. If there is an accidental contact, flush immediately with copious amounts of water and contact a physician if necessary.

With susceptible patients, sensitization to ZAKK[®] Implant cannot be excluded. Do not use ZAKK[®] Implant for patients with allergic reactions against methacrylates. If swallowed accidentally, immediately consult a physician.

Pay attention to the material safety data sheets (MSDS >> www.r-dental.com).

Technical data

Mixing time:	auto-mixing
Working time incl. mixing time (at chairside light, 73 °F):	< 2 min
Light polymerization time (at light sensitive positions):	~ 20 s
Intra-oral polymerization time (self-curing):	> 4 min

ZAKK and KwikkMix are registered trademarks of R-dental Dentalerzeugnisse GmbH.

Caution: For dental use only. US Federal Law restricts this device to sale by or on the order of a dentist or licensed professional, to use or to order the use of this device.

The described product must be used according to the instructions for use and is determined for dental use only. There is no liability for defects resulting from inappropriate use or not intended purposes. The user is responsible for testing the suitability of the product for any purposes not explicitly stated in the instructions for use.

If there are damages to this product within the warranty period, R-dental is obliged to replace or to repair this product. There is no further claim. Outside of the liability by law there is no liability of R-dental for any loss or damages arising from this product, independently whether these damages are direct, indirect, special, incidental or consequential.

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© 2014 R-dental Dentalerzeugnisse GmbH. All rights reserved.
ZAKK y KwikkMix son marcas comerciales registradas de la empresa R-dental Dentalerzeugnisse GmbH.
© 2014 R-dental Dentalerzeugnisse GmbH. Todos los derechos reservados.