VOCO Admira® Protect

CE 0482

EN Instructions for use
MD EU Medical device

Product description:

Light-curing, fluoride-releasing ORMOCER®-based preparation for the treatment of dentine hypersensitivity.

Admira Protect completely reduces hypersensitivity, in most cases after only one application, and additionally protects against external influences and abrasion.

Indications:

- Treatment of hypersensitive dentine
- Treatment of cervical areas (including crown margins, after tooth cleaning or tartar removal, exposed cervical areas following surgery etc.)

Method of use:

1. Preparation

Clean the teeth to be treated.

Isolate areas to be treated with cotton rolls. Remove excess moisture with an oil-free airjet. Do not overdry dentine. The dentine surface should be slightly moist.

Avoid any contamination with blood or saliva after cleaning.

2. Dispensing

When using the bottle:

Dispense **Admira Protect** onto a **Single Tim** or a disposable brush. Alternatively, dispense **Admira Protect** onto a mixing pallet.

When using the *SingleDose*:

Detach a *SingleDose* at the perforated line and turn the printed side upwards. Pierce the foil with a **Single Tim**, enlarge the gap and moisten the brush with liquid through circular motion.

3. Application

Apply **Admira Protect** with **Single Tim** or a disposable brush evenly on all dentine surfaces to be treated and allow to act for 20 s. Then disperse **Admira Protect** with a faint airjet and light-cure with a conventional polymerisation device for 10 s. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. Apply a second layer of **Admira Protect**, disperse it with a faint airjet and light-cure for 10 s.

4. Removal of the inhibition layer

Remove the oxygen-inhibited layer (soft surface) with a cotton pellet or a cotton roll. Check the gingival margins and, if necessary, remove any residual material.

Information, precautionary measures:

Phenolic substances (e.g. eugenol-containing materials) disturb polymerisation and must not be brought into contact with **Admira Protect**. **Admira Protect** contains Bis-GMA, HEMA, BHT, acetone, ORMOCER®s and organic acids. **Admira Protect** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.

Avoid contact with the mucosa. Upon contact with eyes flush immediately with plenty of water and seek medical advice.

Storage:

Close the **Admira Protect** bottle tightly after use. Store at temperatures of 4 °C to 23 °C (storage at 5 °C - 8 °C prolongs the shelflife of the preparation). Do not expose **Admira Protect** to direct sunlight or the chair light.

Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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ORMOCER® is a registered trademark of Fraunhofer Gesellschaft zur Förderung der angewandten Forschung e.V., München, Germany.
VOCO GmbH