VOCO**Bifix® Temp****GB**

Instructions for use

CE 0482**Product description:**

Bifix Temp is a dual-curing, composite-based temporary luting material. The balanced translucency of **Bifix Temp** enables highly aesthetic results to be achieved for non-permanent luting of temporary restorations made from crown and bridge materials as well as temporary CAD/CAM-fabricated pieces. The low film thickness ensures precise placement of temporary restorations at all times. The stability of **Bifix Temp** facilitates particularly simple and reliable removal of excess material.

Working time: > 1,5 min (23 °C)
Setting time: < 4 min (37 °C)

Indications:

- Non-permanent luting of temporary restorations.

Caution:

The surface of composite build-ups and restorations must be coated with a thin layer of insulation (glycerine gel, Vaseline) in order to prevent **Bifix Temp** from luting permanently.

The use of (meth)acrylate-based bonding agents and protective varnishes prior to luting with **Bifix Temp** must be avoided as this can result in a permanent bond.

Application:**1. Preparation:**

Ensure that the surface is clean and dry; if necessary it should be cleaned and dried. Follow the instructions for application on composite build-ups and restorations (see **caution** above).

2. Placing the restoration:

Remove the protective cap from the QuickMix syringe. Due to technical reasons, there may be slight differences in the filling levels of the syringe chambers prior to first use. Therefore remove the syringe cap and extract material until both chambers deliver uniformly. Then fit a mixing tip type 9. The notches on the syringe and the mixing tip must line up exactly. Turn the mixing tip 90° clockwise to lock the tip into place. Apply **Bifix Temp** to the dry temporary restoration and/or the prepared tooth surface. Then place the restoration by applying gentle pressure, generating excess material in the process. The working time at room temperature is at least 1,5 min. Leave the used mixing tip on the syringe to seal it until the next application.

3. Removal of excess material:

Any excess material can be removed using the following methods:

- Wait until the material turns elastic (monitor change in consistency as a consequence of intraoral chemical curing) and remove excess material using a suitable instrument.
- Alternatively, the surplus material can be rendered elastic after short curing with a polymerisation light (around 2 s) and then removed using a suitable instrument.
- If not cured, **Bifix Temp** can be removed using a foam pellet or a similar object.

Removal of excess material following brief curing is recommended. The restoration must be held in position during the polymerisation phase. Contact points and interproximal areas must be cleaned with dental floss. There must be no excess material remaining in the sulcus as this can cause gingival irritation.

4. Curing:

Self-curing is completed 4 min after insertion.

The polymerisation time can be controlled by additional light-curing.

For this, each side is cured for a minimum of 10 s using a polymerisation light with an intensity of at least 500 mW/cm².

Information, precautionary measures:

Bifix Temp is eugenol-free and well-suited for use prior to subsequent adhesive permanent luting as it has no influence on the curing of the materials used.

Only apply **Bifix Temp** using the designated mixing tip, manual mixing is not advisable due to oxygen inhibition.

If, during the course of treatment, the temporary restoration is to be luted again, the application of a thin layer of, for example, glycerine gel or Vaseline to the restoration prior to luting can make it easier to subsequently remove any **Bifix Temp** residue from the temporary restoration. **Bifix Temp** contains dimethacrylates and diacrylates, benzoyl peroxide, amines and BHT. **Bifix Temp** should not be used in case of known hypersensitivities (allergies) to any of these ingredients. In case of contact with the eyes, immediately rinse with plenty of water and seek medical advice.

The warranted product characteristics can only be achieved when using the original Sulzer Mixpac mixing tips as approved by **VOCO**. All the mixing tips for use with our respective products are available from **VOCO**.

Storage:

Store at temperatures of 4 °C to 23 °C.

Do not use after expiry date.

Refrigeration is recommended.

Allow **Bifix Temp** to reach room temperature before use!

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1055 E1 0817 99 © by VOCO
VOCO GmbH