everStick®PERIO

Fibre reinforcement for periodontal splinting

951335 - NU6907
30000531

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gcamerica.com

MADE IN FINLAND

Fibre type: Silanated E-Glass fibre impregnated with Bis-GMA and PMMA

Form: Unidirectional fibre bundle

Diameter: ~ 1.1-1.3 mm

Keep away
from sunlight

Temperature
limit

Caution: US Federal law restricts this device to sale by or on the order of a dentist.

everStick®PERIO

- EN FIBRE REINFORCEMENT FOR PERIODONTAL SPLINTING
- BG ПОДСИЛЕНИ ГЛАС-ФИБРО ВЛАКНА ЗА ПАРОДОНТАЛНО ШИНИРАНЕ
- CS VÝZTUŽ SKLENĚNÝMI VLÁKNY EVERSTICK®PERIO PRO PERIODONTÁLNÍ DLAHOVÁNÍ
- HR STAKLENA VLAKNA ZA POJAČANJE PARODONTNIH UDLAGA
- HU ÜVEGSZÁLAS MEGERŐSÍTÉS A PARODONTÁLIS SÍNEZÉSHEZ
- PL WZMOCNIENIE Z WŁÓKNA SZKLANEGO DO SZYNOWANIA PERIODONTOLOGICZNEGO
- RO FIBRE DE STICLĂ PENTRU CONSOLIDAREA IMOBILIZĂRILOR PERIODONTALE
- RU АРМИРУЮЩЕЕ СТЕКЛОВОЛОКНО ДЛЯ ПАРОДОНТАЛЬНОГО ШИНИРОВАНИЯ
- SK VÝSTUŽ ZO SKLENÝCH VLÁKIEN PRE PARODONTÁLNE DLAHOVANIE
- SL OJAČITVE IZ STEKLENIH VLAKEN ZA PARADONTALNE POVEZAVE
- SR STAKLENA VLAKNA ZA POJAČANJE PARODONTALNIH UDLAGA
- UK СКЛОВОЛОКОННИЙ ЗМІЦНЮЮЧИЙ ЕЛЕМЕНТ EVERSTICK® PERIO ДЛЯ ПЕРІОДОНТАЛЬНОГО ШИНУВАННЯ
- TR FIBER DESTEK PERIODONTAL SPLINTLEME İÇİN
- LV STIKLA ŠĶIEDRU ARMATŪRA PERIODONTĀLAJAI ŠINĒŠANAI
- LT STIKLO PLUOŠTO JUOSTELĖS PERIODONTINIŲ DANTŲ SUTVIRTINIMUI
- ET KLAASFIIBERTUGEVDUS PERIODONTAALSEKS LAHASTAMISEKS

GC

EN	Fibre reinforcement for periodontal splinting	5	SK	Výstuž zo sklených vlákien pre parodontálne dlahovanie	55
BG	подсилени глас-фибро влакна за пародонтално шиниране	11	SL	Ojačitve iz steklenih vlaken za paradontalne povezave	61
CS	Výstuž skleněnými vlákny everStick®PERIO pro periodontální dlahování	18	SR	Staklena vlakna za pojačanje parodontalnih udlaga	66
HR	Staklena vlakna za pojačanje parodontnih udlaga	24	UK	Скловолоконний зміцнюючий елемент everStick® PERIO для періодонтального шинування	72
HU	Üvegszálás megerősítés a parodontális sínézéshez	30	TR	Fiber destek periodontal splintleme için	79
PL	Wzmocnienie z włókna szklanego do szynowania periodontologicznego	36	LV	Stikla šķiedru armatūra periodontāļajai šinēšanai	84
RO	Fibre de sticlă pentru consolidarea imobilizărilor periodontale	42	LT	Stiklo pluošto juostelės periodontinių dantų sutvirtinimui	89
RU	Армирующее стекловолокно для пародонтального шинирования	48	ET	Klaasfiibertugevdus periodontaalseks lahastamiseks	94

Prior to use, carefully read the instructions for use

WHAT IS everStick®PERIO?

everStickPERIO fibre reinforcement is a combination of glass fibres and a permeable polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which facilitates handling of the fibres. The fibre bundle is flexible and sticky, which allows it easily and reliably bond to teeth. The primary application for everStickPERIO fibre reinforcement is the splinting of teeth.

Periodontal surface-retained and intra-coronal splints

- Lingual/palatal splints
- Labial splints
- Occlusal splints

CONTRA-INDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

INSTRUCTIONS FOR USE:

IMPORTANT: The everStickPERIO fibres should be positioned close to the incisal edge to minimise the forces the splint will be subjected to. Also, the splint should not interfere with the occlusal contacts e.g. on the palatal surface of the upper anterior teeth.

Spot fixation with composite does not provide a sufficient bond between the fibre splint and tooth surface. Bond the fibre splint to the teeth for their entire length. Cover the fibre bundle with a thin (0.5 mm)

layer of composite, including the approximal areas, when bonding it to the tooth surfaces. At the occlusal contact in intra-coronal splints, the optimal thickness of the composite layer on top of the fibre is approximately 1–2 mm.

Splints in the upper anterior area do not necessarily have any room palatally for a surface-retained splint, due to the occlusion. Consequently an intra-coronal grooved splint or surface retained labial splint must be considered.

In an anterior area's lingual/palatal splint, a short additional fibre can be used to offset the occlusal forces that tend to loosen fibre splints in the canine region. The additional fibre is attached to the labial surface of the canine and the lateral incisor.

SURFACE-RETAINED PERIODONTAL SPLINTING IN THE ANTERIOR AREA

1. Measure and cut the fibre

Measure the length of fibre needed from the dental arch by using, for example, a periodontal probe or dental floss, to prepare an everStickPERIO fibre splint. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35,6-46,4°F) when you do not need it.

2. Clean the tooth surfaces

The entire length of the fibre splint must be bonded to tooth surfaces. Clean the tooth surfaces with a paste of pumice and water, rinse and air-dry the area. Place wedges in the approximal spaces as necessary, so that the spaces to be cleaned are not filled with composite. If you are working without wedges, be careful not to block these spaces with composite – see item 5.

3. Etch the tooth surfaces

Etch the tooth surfaces and inter-proximal spaces thoroughly with ortho-phosphoric acid, in the area of the splint, in accordance with the instructions of the bonding agent manufacturer. Preferably etch slightly wider than necessary rather than too little. The recommended enamel etching time for surface-retained areas is 45 to 60 seconds.

Rinse with water and air-dry the tooth surfaces thoroughly after etching. As with all bonded restorations a dry operating field is absolutely necessary and rubber dam isolation is highly recommended.

4. Bond the tooth surfaces

Use the adhesive bonding technique for bonding teeth according to the instructions of the bonding agent manufacturer. Apply the bonding agent to the entire area to be bonded. Light-cure the bonding agent as described by the manufacturer.

5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-ænial Universal Flo) on to the surface of the teeth for the entire width of the fibre bundle. Carefully cover the bonding area with a thin layer (about 0.5 mm)

of composite including the approximal spaces. Leave enough space for cleaning the approximal spaces. Do not cure the composite during this phase.

6. Position and light-cure the fibre

Remove the white protective paper and use tweezers to pick the fibre up from the silicone groove. Remove any residual silicone from the fibre bundle. Position the fibre bundle on top of the uncured flowable composite. Aim to place the fibre as incisally as possible in the anterior area. Make sure that it will not be in occlusion. Position one end of the fibre bundle first by pressing it down with the StickSTEPPER instrument (sterilize the instrument before use). Pre-cure the fibre in place, one tooth at a time, for about five seconds, using a curing light. The wide-tipped StickSTEP-

PER instrument shields the rest of the fibre from light. As the fibre is made of light-conducting material, it is recommended to direct the light-curer away from the uncured fibre bundle. Press the fibre into the approximal spaces as well. Make sure that the embrasures are not blocked with fibre and composite.

7. Cover and finish the splint

After pre-curing, cover the entire fibre splint with a thin layer of composite. Note that in surface-retained areas the fibre bundle can be coated with a thin layer (0.5 mm) of composite. Then light cure the whole splint for 40 seconds, one tooth or coverage area of the light-curer at a time. Be careful not to cut the fibre when finishing/polishing the splint.

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The steps for the intra-coronal splint are the same as for the surface retained splint except for the preparation of a groove, mesial to distal, in the teeth to be splinted.

1. Prepare the groove for the teeth

Prepare a groove for the teeth to be splinted with a minimum width of 2 mm. Ideally, it should stay within the enamel as this will provide the best bond. At the occlusal contact, the optimal thickness of the composite layer on top of the fibre is 1–2 mm. The margins of the groove should be beveled as this will ensure the best marginal integrity and enlarge the etched enamel surface for composite bonding. Place wedges in the approximal spaces, if possible, so that the embrasures are not

filled with composite. If you are working without wedges, note that the embrasures must remain free of composite.

2. Measure and cut the fibre

Measure the length of fibre needed to prepare an everStickPERIO fibre splint for the prepared groove by using, for example, a periodontal probe or dental floss. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35,6-46,4°F) when you do not need it.

3. Etch the teeth

Etch the prepared groove thoroughly with ortho-phosphoric acid in accordance with the instructions of the bonding agent manufacturer. Rinse with water and air-dry the tooth surfaces thoroughly after etching. As with all bonded restorations a dry operating field is absolutely necessary and rubber dam isolation is highly recommended.

4. Bond the teeth

Bond the prepared groove area in accordance with the instructions of the bonding agent manufacturer. Light-cure the bonding agent as described by the manufacturer.

5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-aenial Universal Flo) into the prepared groove. Be careful not to block the

embrasures with composite. Do not light cure the composite during this phase.

6. Position and light cure the fibre

Remove the white protective paper and use tweezers to pick the fibre up from the silicone groove. Remove any residual silicone from the fibre bundle. Place the fibre bundle into the groove on top of the uncured flowable composite.

Press the fibre bundle into the flowable composite with the StickCARRIER instrument. Pre-cure the fibre in place, one tooth at a time, for about five seconds, with a curing light. During light curing use the wide-tipped StickSTEPPER instrument to shield the rest of the fibre from light. As the fibre is made of light-conducting material, it is recommended to direct the light-curer

away from the uncured fibre bundle.

7. Cover and finish the splint

Cover the splint and fill the groove with a 1–2 mm composite layer. Then light cure the whole splint for 40 seconds, one tooth or coverage area of the light-curer at a time. Adjust the occlusion and finish the splint. Be careful not to cut the fibre when finishing/polishing the splint.

TIPS AND RECOMMENDATIONS

- The fibre must always be covered entirely with composite.
- Use rubber dam isolation to ensure a dry operating field.
- Always try to handle the fibres with instruments to avoid contamination e.g. powdered gloves.
- Use the StickSTEPPER instrument for

sectional curing of the fibre bundle.

- Position the splint as close as possible to the incisal edges of the teeth in the anterior area.
- If, after placing the fibre, you notice it is too long, shorten it with a diamond bur during the finishing phase of the splint. Apply some enamel resin to the exposed fibre surface (e.g., StickRESIN), blow the resin into a thin layer and light cure. Cover the fibre again carefully with composite.
- Avoid cutting the fibres during the finishing phase.
- At the occlusal contact area, the optimal thickness of the composite layer on top of the fibre splint is approximately 1–2 mm.

STORING: everStick products should always be stored in a refrigerator (2-8°C/35,6-46,4°F). In addition, the products should be protected

from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

(Shelf life : 2 years from date of manufacture)

PACKAGES

everStick Intro

8cm everStickC&B; 8cm everStickPERIO;

1x StickREFIX D silicone instrument
everStick COMBI
 8cm everStickC&B; 8cm everStickPERIO;
 30 cm² everStickNET; 5x everStickPOST Ø 1.2;
 1x StickSTEPPER hand instrument;
 1x StickREFIX D silicone instrument

everStick Starter Kit:

8 cm everStickC&B; 8 cm everStickPERIO;
 30 cm² everStickNET; 5 x everStickPOST
 Ø 1.2; 5 ml StickRESIN bottle; 2 ml G-ærial
 Universal Flo syringe; 20 dispensing tips and
 1 light protective cap; 1x StickREFIX D
 silicone instrument; 1x StickSTEPPER;
 1x StickCARRIER

Refills

2 x 12 cm fibre bundle
 1 x 8 cm fibre bundle

NOTE: Products should be used clinically with care and the patient should be warned not to abrade the fitting surface so as to avoid exposing

irritation-causing fibers.

If the surface of the fibre bundle feels dry, but it is fully bendable and not polymerized, adding a drop of resin (such as StickRESIN) will return the flexibility/workability of the material.

Polymerization can be observed as white spots at bending area when bending the bundle. The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will still continue during the next 24 hours. StickSTEPPER, StickCARRIER hand instruments and StickREFIX D, StickREFIX L silicone instruments must be sterilized before use.

WARNING: Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes

in contact with resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerized everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases. The use of powder free gloves is recommended with everStick products. Polymerize everStick before waste disposal.

US Federal law restricts this device to sale by or on the order of a dentist.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<http://www.gceurope.com> or for The Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

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CEP: 02022-051, CNPJ: 08.279.999/0001-61, RESP. TÊC: MAYARA DE SANTIS

RIBEIRO - CRO/SP 105.982

RESPONSIBLE MANUFACTURER IN CANADA

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Prior to use, carefully read the instructions for use.

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SURFACE-RETAINED PERIODONTAL SPLINTING IN THE ANTERIOR AREA

1. Measure and cut the fibre

Measure the length of fibre needed from the dental arch by using, for example, a periodontal probe or dental floss, to prepare an everStickPERIO fibre splint. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35.6-46.4F) when you do not need it.

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5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-aenial Universal Flo) on to the surface of the teeth for the entire width of the fibre bundle. Carefully cover the bonding area with a thin layer (about

0.5 mm) of composite including the approximal spaces. Leave enough space for cleaning the approximal spaces. Do not cure the composite during this phase.

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1. Prepare the groove for the teeth

Prepare a groove for the teeth to be splinted with a minimum width of 2 mm. Ideally, it should stay within the enamel as this will provide the best bond. At the occlusal contact, the optimal thickness of the composite layer on top of the fibre is 1–2 mm. The margins of the groove should be beveled as this will ensure the best marginal integrity and enlarge the etched enamel surface for composite bonding. Place wedges in the approximal spaces, if

possible, so that the embrasures are not filled with composite. If you are working without wedges, note that the embrasures must remain free of composite.

2. Measure and cut the fibre

Measure the length of fibre needed to prepare an everStickPERIO fibre splint for the prepared groove by using, for example, a periodontal probe or dental floss. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35.6-46.4F) when you do not need it.

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5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-aenial Universal Flo) into the

prepared groove. Be careful not to block the embrasures with composite. Do not light cure the composite during this phase.

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Press the fibre bundle into the flowable composite with the StickCARRIER instrument. Pre-cure the fibre in place, one tooth at a time, for about five seconds, with a curing light. During light curing use the wide-tipped StickSTEPPER instrument to shield the rest of the fibre from light. As the fibre is made of light-conducting material,

it is recommended to direct the light-curer away from the uncured fibre bundle.

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TIPS AND RECOMMENDATIONS

- The fibre must always be covered entirely with composite.
- Use rubber dam isolation to ensure a dry operating field.
- Always try to handle the fibres with instruments to avoid contamination e.g. powdered gloves.

- Use the StickSTEPPER instrument for sectional curing of the fibre bundle.
- Position the splint as close as possible to the incisal edges of the teeth in the anterior area.
- If, after placing the fibre, you notice it is too long, shorten it with a diamond bur during the finishing phase of the splint. Apply some enamel resin to the exposed fibre surface (e.g., StickRESIN), blow the resin into a thin layer and light cure. Cover the fibre again carefully with composite.
- Avoid cutting the fibres during the finishing phase.
- At the occlusal contact area, the optimal thickness of the composite layer on top of the fibre splint is approximately 1–2 mm.

STORING: everStick products should always be stored in a refrigerator (2-8°C/35.6-46.4F). In addition, the products should be protected from light by packing them in the sealed

foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

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Refills

2 x 12 cm fibre bundle
1 x 8 cm fibre bundle

NOTE: Products should be used clinically
with care and the patient should be warned
not to abrade the fitting surface so as to
avoid exposing irritation-causing fibers.
If the surface of the fibre bundle feels dry, but it
is fully bendable and not polymerized, adding

a drop of resin (such as StickRESIN) will return
the flexibility/workability of the material.
Polymerization can be observed as white spots
at bending area when bending the bundle.
The everStick fibres do not achieve their full
strength immediately after the final light-curing
of 40 seconds. The polymerization of the fibres
will still continue during the next 24 hours.
StickSTEPPER, StickCARRIER hand instru-
ments and StickREFIX D, StickREFIX L silicone
instruments must be sterilized before use.

WARNING: Personal protective equipment
(PPE) such as gloves, face masks and safety
eyewear should always be worn. Unpolymer-
ised resin can cause skin sensitisation to
acrylates in some people. If your skin comes
in contact with resin, wash it thoroughly with
soap and water. Avoid contact of uncured
material with skin, mucous membrane, or eyes.
Unpolymerized everStick products may have

a slight irritating effect and lead to sensitiza-
tion to methacrylates in rare cases. The use of
powder free gloves is recommended with
everStick products. Polymerize everStick
before waste disposal.

US Federal law restricts this device to sale by
or on the order of a dentist.

Some products referenced in the present IFU
may be classified as hazardous according to GHS.
Always familiarize yourself with the safety
data sheets available at:

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