# exaplast, exasoft, exagel

**GB**  
Instructions for use

Precision impression material, silicone based, condensation curing

Fig. 1

## 1. Tray preparation

Select a tray type suited for desired application. For optimal adhesion of impression material, we recommend use of **sili Adhesive liquid** or **spray**. Apply a thin layer of adhesive film onto the tray and let dry for 60 seconds (Fig. 1). Take preliminary impression with **exaplast putty**.

Fig. 2

## 2. Preliminary impression

Fill measuring scoop evenly to top with **exaplast putty** (Fig. 2). Smooth **exaplast putty** onto mixing pad with a spatula and make a series of crosswise grooves (Fig. 3). Add recommended amount of **exagel cat** (Fig. 4). Mix thoroughly, first using spatula and then by kneading between fingers until a uniform colouring is achieved (approx. 30 sec.). To facilitate the flowing out of excessive light bodied material from the impression, cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended. Clean the surface of the preliminary impression prior to the corrective impression.


Fig. 3

## 3. Corrective impression

Extrude **exasoft wash** and **exagel cat** in equal lengths and thickness onto mixing pad (Fig. 5). Mix both components for approx. 30 sec. with spatula until a homogeneous uniform colouring is achieved. Load preliminary impression (Fig. 6). A normal impression syringe can be used for bubble free application of **exasoft wash** (Fig. 7).

Fig. 4

## 4. Disinfecting

After removing tray from mouth, rinse impression under lukewarm water. Then impression can be disinfected by immersing into 2% glutaraldehyde for 15 min.


Fig. 5

## 5. Model casting

Casting between 30 minutes and 72 hours after removal of the impression. All class III and IV dental plasters and standard model-casting acrylics can be used.


Fig. 6

## 6. Electro plating

The impression can be electro plated using all commercially available copper or silver plating baths.

Fig. 7

## 7. Cleaning the tray


After setting, material can be dislodged with a blunt instrument. Remaining adhesive film can easily be removed with **sili Solvent** spray (FCHC-free; use only in well ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

## Important working hints


- We recommend using **exasoft wash** for corrective impression after making a preliminary impression with **exaplast putty**.
- Overdosing of Activator and increased temperatures accelerate, underdosing and decreased temperatures retard the setting time.
- Activator tubes should be closed tightly immediately after use.
- **exasoft wash**: Eventually separated silicone oil during storage of the product also can be mixed, without having any negative effect on the working process.
- The setting reaction is not influenced by that.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Avoid contact of Activator with eyes, may cause irritation. If necessary, flush eyes with plenty of water and seek medical advice immediately.
- Don't leave any residual material in the patient's mouth.
- Please follow the instructions of the safety data sheets!

Avoid contact with eyes and skin!  
Wear suitable eye protection and gloves!

## exaplast putty / exagel cat

## exasoft wash / exagel cat

## Indications for use:

- **exaplast putty**  
Preliminary impression material for two-stage impressions, putty for double mix technique, situation impressions, check bite
- **exasoft wash**  
Corrective material for double impression and double mix technique, relining impressions

## Technical data:

**exaplast putty**  
DIN EN ISO 4823 – Type 0

**exasoft wash**  
DIN EN ISO 4823 – Type 3

- **Dosage:**  
**exaplast putty/ exagel cat**  
1 scoop = 13.5 g /  
6 cm = 0.3 g  
**exasoft wash / exagel cat**  
10 cm = 5.0 g /  
10 cm = 0.5 g
- **Colour code:**  
exaplast putty: yellow  
exasoft wash: pink  
exagel cat: blue
- **Mixing time:**  
approx. 30 sec., each
- **Working time:**  
approx. 1 min. 15 sec. \*, each
- **Setting time:**  
exaplast putty:  
approx. 4 min. 45 sec. \*  
exasoft wash:  
approx. 3 min. 45 sec. \*
- **Time in the mouth:**  
exaplast putty:  
approx. 3 min. 30 sec.  
exasoft wash:  
approx. 2 min. 30 sec.
- **Strain in compression:**  
exaplast putty: 2.9 %  
exasoft wash: 7.0 %
- **Recovery from deformation:**  
exaplast putty: 98.0 %  
exasoft wash: 98.5 %
- **Linear dimensional change:**  
exaplast putty: 0.3 %  
exasoft wash: 0.7 %
- **Application:**  
At 23 °C ± 2 °C /  
73 °F ± 4 °F, 50 ± 5 %  
relative humidity

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

## Storage/ shelf life

The materials fulfill the intended purpose at least until the imprinted expiry date in well sealed containers at temperatures of 15-25 °C/59-77 °F and 50 % relative humidity. Impressions should be stored at normal room temperature. Avoid exposure to heat and sun.

## Ordering information:

<b>exaplast</b> Set Putty, jar of 910 ml Wash, tube of 140 ml Cat, tube of 60 ml	<b>03133</b>
<b>exaplast putty</b> Standard packing jar of 910 ml	<b>02907</b>
<b>exasoft wash</b> Standard packing tube of 140 ml	<b>02910</b>
<b>exagel cat</b> tube of 60 ml	<b>03125</b>

Made in Germany

09/2017

**DETX**

GmbH & Co. KG  
Postfach 100225 · 76256 Ettlingen / Germany  
Telefon 07243/510-0 · Fax 07243/510-100  
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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