FREEFORM®

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Lichthärtender, transparenter Modellierkunststoff, MMA-frei, biokompatibel, röntgentransluzent, sterilisierbar


Abb. 1

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
Fig. 6

1. Preparation

Block out plaster model, wet and carefully apply alginate isolation (repeat if required) (Fig. 1). Allow to dry thoroughly. The use of petrolatum is recommended to isolate the counter-bite in the articulator. To avoid sticking, apply a thin coat of petrolatum to fingers, gloves (i.e. nitrile), modelling instruments etc. This considerably facilitates preforming and modelling.

2. Application

Remove an appropriate amount of **FREEFORM® plast** from the tin, reclose tightly, and shape the material to a homogenous ball in your palms applying slight pressure (Fig. 2). Assess the amount of material, homogenise and roll again after adding or removing material to avoid inclusions. Further shaping depends on the prevailing situation (i.e. denture model U-shape or basic ball). **Apply FREEFORM® plast** to the corresponding areas and adapt gently (Fig. 3). Form the desired shape by gentle, repeated pressing with a finger tip. Observe appropriate dimensioning (recommended minimum thickness 1.2 mm). If the antagonists in the articulator are to be included, then apply petrolatum both to the model as well as the **FREEFORM® plast** surface to avoid sticking and deformation. In accessible areas it is possible to smooth the surface by smoothing with finger tips.

3. Curing

FREEFORM® plast can be polymerised in all commercially available light curing devices with a light spectrum of 320 nm - 400 nm. Units with cooling of the polymerisation chamber are also suitable. The average polymerisation time in light-curing devices with UV/UV-A lamps is approx. 10 min. Wet again for 2-3 before removal. This softens the insulation and facilitates removal. Perform loosening and removal of the model with great care and in direction of insertion (risk of breaking if wedged).

Light-cure the opposite side of the model for approx. 5 minutes if subsequent mechanical polishing is intended. If **FREEFORM® coat** is used, this step is not required as a further polymerisation process is required. Please note that the polymerisation time always depends on the lamps in the device. Uniform curing times can only be ensured with regularly serviced devices. Observe the device manufacturer's instructions for use.

4. Processing

Following complete polymerisation of **FREEFORM® plast**, remove the slightly sticky inhibition layer by rubbing with isopropanol.

If auxiliary parts, such as implant drill sleeves or X-ray reference spheres are to be integrated, then these can be worked in with **FREEFORM® fixgel** after placing (Fig. 4). To this purpose, roughen the surface in the surrounding area, apply **FREEFORM® fixgel** and polymerise for 5-10 min. If this is to be performed on already completed, older templates, then it is essential to apply **FREEFORM® bond** first, followed by **FREEFORM® fixgel**, and then cure both together (Fig. 5).

This also allows complementing, reinforcing and repairs to **FREEFORM® plast** and **FREEFORM® fixgel**.

In order to join **FREEFORM® plast** with a deep-drawn splint, ensure that this is suitable for connecting to PMMA resins (i.e. ERKODUR, BIOCRYL). Then roughen the splint by sandblasting, apply **FREEFORM® bond** and allow to work for 1-2 min. Then adapt **FREEFORM® plast** and cure for approx. 10 min. Tungsten carbide burs are suited for processing and grinding. Sand the surface with emery paper and then polish mechanically. **FREEFORM® coat** gloss lacquer can be used in place of manual high gloss polishing (Fig. 6). First additionally smooth the surface by sandblasting (aluminium oxide 100).

Sealing of the top and bottom sides can be performed in two work steps. Tip: only apply the lacquer very thin to the bottom side, possibly disperse with air, so as not to impair the fit. Light-cure **FREEFORM® coat** within approx. 5 min. per side.

5. Disinfection and sterilisation

Disinfection

FREEFORM® can be disinfected with standard agents. Please observe the manufacturer's compatibility data.

Sterilisation

Locally applicable regulations and hygiene guidelines are to be observed.

FREEFORM® is validated for the "W&H Lisa 522, Programme UNIVERSAL 134" (see technical specifications). Steam sterilisation may only be performed with equipment that complies with the standards EN 13060 and EN 285. Validation of the sterilisation processes was conducted according to EN ISO 17664. Responsibility for sterility lies with the user.

Important notes on processing

- **FREEFORM® plast** can be combined with PMMA resins and deep-drawn splints (observe suitability according to manufacturer's data) for stabilising or adapting the surface.
- Observe material thicknesses, minimum 1.2 mm, maximum 5 mm, up to 10 mm for polymerisation on both sides.
- Biocompatibility is only guaranteed with complete polymerisation.
- Avoid high temperatures during polymerisation as this can lead to distortion of the objects. Non-cooled stroboscopic light units are therefore not suitable.
- Prior to polymerisation, check the contoured work for possible lifting or separating from the model.
- Close tin tightly immediately after removal of material.
- Milling and polishing should only be performed on work stations with suction devices.
- Avoid contact with mucous membranes and eyes. In case of accidental contact, rinse with adequate running water, consulting an ophthalmologist if necessary.
- Non-cured material can lead to irritation and sensitisation upon contact with skin or mucous membranes. Avoid contact with skin.
- Wear gloves and eye protection when processing material which has not yet set.
- **FREEFORM®** can be disinfected using impression disinfection agents, for example.
- The product was developed for use in the dental field and may only be used by professionals for its intended use.
- Please follow the instructions of the safety data sheets!

Fields of application:

- Fabrication of individual parts, for example, drilling templates for implant dentistry, X-ray templates, base plates, splints, fixation or transfer keys.

Contra-indication:

Sensitisation through **FREEFORM®** is extremely rare. If such reactions occur, the application of **FREEFORM®** is to be discontinued. **FREEFORM®** may not be dispensed intraorally in non-polymerised condition.

Technical specifications:

Volume/weight

FREEFORM® plast:

60 g tin

FREEFORM® coat:

10 ml bottle

FREEFORM® fixgel:

3 g syringe

FREEFORM® bond:

5 ml bottle

Dispensing & application

Direct application, manual contouring

Product colour

transparent

Curing

Light spectrum of equipment 320 nm - 400 nm

UV/UV-A lamps approx.

10 min. on model + approx.

5 min. bottom side without model

Sterilisation process

W&H Lisa 522, Programme UNIVERSAL 134 with the following parameters:

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