FREEPRINT® cast

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FREEPRINT® cast

Light-curing resin based on (meth)acrylate, for cast parts in dental precision casting, burn-out for DLP printers with UV-LED 385 nm / LED with 405 nm

Suitable for the following DLP-Printer

Asiga 385 nm		405 nm	MiiCraft 385 nm	Rapidshape 385 nm		405 nm	W2P 385 nm
MAX Pro2 Pico2	MAX 405 Pro2 405 Pico2 405		MiiCraft y-Serie	DII-Serie	D-Serie		Solflex Serie

Printers may only be operated using material parameters authorised by DETAX.

Important notes

Only for dental use by qualified personnel.

Processing

- ▶ Phosphate-bonded investment material can be used for fabrication of the casting mould. A pre-heating temperature of approx. 800 °C, is to be maintained for approx. 30-45 minutes, depending on the size of the casting mould, regardless of the final casting temperature.
- ▶ Homogenize/roll the material prior to processing, i.e. with a roller mixer.
- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the object are completely cured (observe process description on page 2).
- ▶ Processing temperature 23 °C ± 2 °C.

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.

Storage

- ▶ **FREEPRINT® cast** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® cast** may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

Adverse effects

Product may cause allergic reactions.

Indication:

Fabrication of residue-free combustible cast objects for precision casting technology employing the DLP process.

Processing:

at 23 °C ± 2 °C

Storage

Ordering information:

FREEPRINT® cast 385
500 g bottle, red **03710**
1000g bottle, red **03789**

FREEPRINT® cast 405
500 g bottle, red **02890**
1000g bottle, red **03788**

FREEPRINT® model 385
1.000 g bottle
sand **03778**
ivory **03780**
grey **03782**

FREEPRINT® gingiva 385
500 g bottle, gingiva **02820**
1.000 g bottle, gingiva **02843**

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 07243/510-0 · Fax: 07243/510-100
www.detax.de · post@detax.de

Made in
Germany
04/2019

More information at
www.detax.de

FREEPRINT® cast

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

Post-curing

Post-curing is performed with a xenon photoflash unit (e.g. Otoflash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

FREEPRINT® cast was developed for use in the dental field and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.

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