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Light-curing resin based on (meth)acrylate, biocompatible, for the generative fabrication of occlusal splints, drilling templates, X-ray templates for DLP printers with UV-LED 385 nm / LED with 405 nm

Suitable for the following DLP-Printer

Asiga 385 nm		405 nm	MiiCraft 385 nm	Rapidshape 385 nm		405 nm	W2P 385 nm		
MAX UV	Pro2 UV	Pico2 UV	MAX 405	Pro2 405	Pico2 405	MiiCraft y-Serie	DII-Serie	D-Serie	Solflex Serie

Printers may only be operated using material parameters authorised by DETAX.

Important notes

This is a medical device, only to be used by trained specialist personnel.

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 2).
- ▶ Homogenize/roll the material prior to processing, i.e. with a roller mixer.
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ Maximum curing depth* at direct post-exposure:
clear: 3,5 mm
*In case of large objects and exposure on both sides, the material thickness can be up to 7 mm (Example FREEPRINT® splint – with a curing depth of 3,5 mm). Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ▶ Seal surface with FREEFORM® coat.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.

Storage

- ▶ FREEPRINT® splint is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

Contraindication

Contains (meth)acrylics and phosphine oxide.

Some ingredients of FREEPRINT® splint may cause allergic reactions in predisposed persons. In such cases refrain from using the product. FREEPRINT® splint only insert intraorally in completely polymerised state.

Adverse effects

Product may cause allergic reactions.


Indication:

Drilling and X-ray templates for implant dentistry, occlusal splints, fixation and transfer keys

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Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath.

Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otofash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

Surface processing

Varnish the surface with FREEFORM® coat

FREEPRINT® splint was developed for use in the dental field and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.