# glassbite

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# glassbite

VPS bite registration material, transparent,  
fast setting – high final hardness


Fig. 1

## 1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Insert mixing tip and rotate clockwise to a stop (Fig. 3). Store cartridge with mixing cannula attached for subsequent use.


Fig. 2

## 2. Instruction for use

Extrude material from cartridge with mixing cannula or with attached contouring tip onto teeth (Fig. 4). At terminal occlusion application from vestibular is possible. Let registration harden in terminal occlusion and remove from the mouth after 1 min.


Fig. 3

## 3. Disinfection

After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde or other usual solutions.

## Important working hints

- Remove any retraction astringents or mouth wash solutions prior to the registration – could affect setting reaction.
- Do not intermix with condensation curing silicones.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Remove all residues of registration material from the patient's mouth.
- Please follow the instructions of the safety data sheet!


Fig. 4

## Further information:

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

## Indications for use:

- Bite registration with crown and bridge impressions and for orthodontic and gnathologic indications, key for intraoral registration
- Application in the dental office and the laboratory as a transparent matrix material for light-curing composites
- For fixation of X-ray reference spheres, i.e. for implant planning with digital 3D X-ray or digital volume tomography

## Technical Data:

- Mixed volume:** 50 ml (cartridge)
- Mixing ratio:** 1:1
- Colour code:** base: transparent catalyst: transparent
- Mixing time:** omitted (Automix2-system)
- Working time:** approx. 30 sec.\*
- Setting time:** approx. 1 min. 30 sec.\*
- Time in mouth:** approx. 1 min.
- Hardness:** approx. 80 Shore A
- Linear dimensional change:** approx. 0.2 %
- Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

## Ordering information:

- glassbite Standard packing 02571**  
1 x 50 ml cartridge  
4 mixing cannulas, green  
4 contouring tips for mixing cannulas green
- Multipack 02572**  
4 x 50 ml cartridges  
16 mixing cannulas, green  
16 contouring tips for mixing cannulas green
- Mixing cannulas 02705**  
green, 48 pcs.
- Contouring tips 02648**  
for mixing cannulas green  
12 pcs.

**DETAX**

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 072 43/510-0 · Fax: 072 43/510-100  
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personell).

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