iCEM® Self Adhesive

- | | |
|--|---|
| (DE) Gebrauchsanweisung (GB) Instructions for use (FR) Mode d'emploi (ES) Instrucciones de uso (IT) Istruzioni per l'uso (PT) Instruções de uso (NL) Gebruiksaanwijzing (SE) Bruksanvisning (DK) Brugervejledning | (NO) Bruksanvisning (FI) Käyttöohjeet (GR) Οδηγίες χρήσης (CZ) Návod k použití (HU) Használati utasítás (LV) Lietošanas instrukcija (LT) Naudojimo instrukcija (PL) Instrukcja obsługi |
|--|---|

KULZER MITSUI CHEMICALS GROUP

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Made in USA

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iCEM® Self Adhesive

Instructions for use (GB)

ICEM Self Adhesive
Dual curing self-adhesive composite cement, radio-opaque

ICEM Self Adhesive is a dual curing, radio-opaque, fluoride releasing, self-adhesive composite cement for the cementation of indirect restorations.

Compositzement
iCEM Self Adhesive is based on acidified urethane and di-, tri-, and multifunctional acrylate monomers and contains 49% filler by weight of sub-micron and micron sized particles.

Clinical Indication

• Luting metallo-ceramic and metal (both precious and non-precious) crowns and bridges

• Luting crowns and bridges using reinforce ceramics (e.g. zirconia)

• metal, ceramic and composite onlays

• precious and non-precious metal, titanium and fiber posts

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Clinical treatment protocol

Ceramic, metal, and cured composite surfaces should be dry.

All-Ceramic Restoration

Follow ceramic restoration manufacturer's recommended tooth-relocation specifications when preparing teeth for reinforced all-ceramic restorations. Reinforced all-ceramic restorations should be closely fitted. Restrict the laboratory to use of space to accommodate cement with a low film thickness. To optimize mechanical retention, a well-etched or roughened ceramic surface is recommended. Optimal mechanical retention is indicated for bonding to zirconia type surfaces.

Instructions for Crown and Bridge, Inlay and Onlays

1. The dentin should etch or abrade the tooth contacting surfaces of restorations to be cemented. If not, etch or microabrade these surfaces. Rinse and dry.

2. Rinse tooth and dry. Remove excess water. LEIHZAHN SURFACES SLIGHTLY MOIST.

3. Etching and bonding agents are required not dentin and cut enamel. Etching uncut enamel is indicated.

4. If using a ceramic, metal and composite in the same preparation, clean, etch or mechanically roughen the surface, rinse and dry.

5. Isolate adjacent teeth.

6. Remove cap from the double barrel syringe. If necessary bleed the syringe so that base and catalyst are at the orifice of the syringe barrels. Place a mixing tip on the automic syringe.

7. To ensure homogeneous mixing of the two components, tighten the screw on the first 2-3 mm of the dispersion material only during the first use of the automic syringe.

8. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

9. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

10. Maintain positive pressure on the restoration for 2.5 minutes. During this period, the excess material at the margin of the restoration is removed with a brush or with a suitable instrument after brief light polymerization for 1-2 seconds. Contact and interproximal areas with floss. Cement may bond to adjacent teeth if excess is not removed.

11. Maintain positive pressure on the restoration for 2.5 minutes. During this period, the excess material at the margin of the restoration is removed with a brush or with a suitable instrument after brief light polymerization for 1-2 seconds. Contact and interproximal areas with floss. Cement may bond to adjacent teeth if excess is not removed.

12. Verify no excess cement remains on the restoration during 2 minutes. Periodically remove the excess of material from the margins of the restoration with a small amount of water. Cement may bond to adjacent teeth if excess is not removed.

13. Remove the excess cement from the restoration during 2 minutes. Periodically remove the excess of material from the margins of the restoration with a small amount of water. Cement may bond to adjacent teeth if excess is not removed.

14. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

15. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

16. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

17. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

18. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

19. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

20. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

21. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

22. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

23. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

24. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

25. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

26. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

27. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

28. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

29. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

30. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

31. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

32. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

33. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

34. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

35. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

36. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

37. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

38. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

39. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

40. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

41. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

42. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

43. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

44. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

45. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

46. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

47. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

48. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

49. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

50. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

51. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

52. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

53. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

54. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

55. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

56. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

57. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

58. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

59. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

60. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

61. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

62. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

63. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

64. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

65. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

66. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

67. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

68. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

69. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

70. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

71. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

iCem® Self Adhesive
Densitivende, selvadherende, rentgenkopositive kompositcement

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