# implantlink® semi Classic

DE

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# implantlink® semi Classic

**GB**  
Instructions for use

## Resin-based temporary luting cement for implant-retained restorations

Dual-curing semi-permanent implant cement with very low displacement resistance and thus extremely low film thickness. Due to highly cross-linked, ductile polymer structure, high compression strength and hermetic sealing are achieved. Excess material can easily be removed in large pieces (no crumbling) during gel phase. Luted restorations are effortlessly removed without damage with standard devices (Crown Butler, Corona Flex). Compatible with all conventional material combinations, i.e. zirconium and metal frameworks, ceramics, titanium. Antibacterial, eugenol-free, odourless and tasteless; natural opaque colour for screening metallic abutments.


Fig. 1

### 1. Preparation and cleaning of implant abutments and inner surfaces of restorations

Carefully check fit, occlusion and proximal contact points. Complete all finishing procedures. Degrease, clean (e.g. alcoholic solutions) and dry prepared abutments. Make sure the interior surfaces of the superstructure are clean and dry. (Fig. 1).


Fig. 2

### 2. Mixing and applying

Squeezing ensued by mini-mix system. Insert the plunger in the body of the cartridge. Remove the cap from the cartridge by turning anti-clockwise. Bleed the syringe before installation of the automix cannula until material emerges uniformly from both openings (Fig. 2). Mount the mixing cannula by using the guides on the cannula and the cartridge and lock in place by rotating in the opposite direction (Fig. 3). Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed homogeneously (Fig. 4). Only then dosing can be performed individually. Leave mixing cannula on cartridge after use, thus sealing the cartridge.


Fig. 3

### 3. Inserting the restoration

Apply a thin layer of **implantlink® semi Classic** directly into the superstructure (Fig. 5). Carefully position the restoration on the implant abutment with slight pressure (Fig. 6); then press down firmly. Remove excess material during the gel phase, after around 2 - 3 minutes, with a suitable instrument (Fig. 7).


Fig. 4

### 4. Working time / light curing

The working time at room temperature is around 60 seconds. The cement is cured at body temperature in about 5 - 6 minutes. Direct exposure to light accelerates setting time. Light curing (marginal region, excess material) is achieved with approx. 20 seconds per surface (Fig. 8). When **implantlink® semi Classic** is used with translucent materials (e.g. thin layers of ceramic provisional crowns and bridges), increase exposure to light through the material to 30 - 40 seconds per surface.

Fig. 5

### Important working hints

- Do not rough or sand-blast abutment surface.
- Do not damage the underlying titanium surface when removing excess material. If necessary, use plastic instruments (implant scalers).
- Leave no cement residues in the mouth or gingival sulcus.
- Use no adhesion-reducing additives. They could uncontrollably render the retentive qualities of **implantlink® semi Classic**.
- Do not manually remove or mix cement components.
- Avoid contact of unmixed components with skin or mucosa. Wear suitable eye protection and gloves!
- Do not use **implantlink® semi Classic** for the luting of hybrid restorations with natural teeth and implant abutments.
- Mounting should be checked after approx. 3-6 months. Then insert the restoration definitively or temporarily again.
- Please follow the instructions on the safety data sheet!


Fig. 6

Only for dental use by qualified personnel.



Fig. 7

### Autopolymerisation:

### Light-curing:


Fig. 8

# DETAX

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100  
www.detax.de · post@detax.de

more information  
<https://youtu.be/pyKtNdyfa54>

Made in Germany

12/2017

### Indication for Use:

- Temporary luting of implant-retained restorations

### Contra-indication:

- In very rare cases **implantlink® semi Classic** may cause sensitization. If such reactions should occur, discontinue using **implantlink® semi Classic**.

### Technical data:

- **Mixed volume:** 5 ml (mini-mix)
- **Mixing ratio:** 4:1
- **Colour code:** Base: white-opaque Catalyst: semi-transparent
- **Mixing time:** None, (mini-mix system)
- **Working time:** approx. 60 sec.
- **Gel phase:** approx. 2 - 3 min.\* (for removal of excess material during placement)
- **Setting time in the mouth:** approx. 5 - 6 min.\*
- **Setting time with light-curing:** approx. 20 sec. per surface; 30-40 sec. through translucent materials
- **Film thickness:** approx. 10 µm
- **Compression strength:** approx. 85 MPa
- **Application temperature:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:** 2 °C / 36 °F, 25 °C / 77 °F

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Higher temperatures accelerate, lower temperatures retard a. m. times.

### Ordering information:

**implantlink® semi Classic**  
Standard packing 03092  
cartridge of 5 ml  
mini-mix 4:1  
10 mixing cannulas, brown

**implantlink® semi Forte**  
Standard packing 03371  
cartridge of 5 ml  
mini-mix 4:1  
10 mixing cannulas, brown

**implantlink® semi Xray**  
Standard packing 02195  
cartridge of 5 ml  
mini-mix 4:1  
10 mixing cannulas, brown

Mixing cannulas 4:1 02591  
brown, 25 pcs.

**monoprint supra**  
Transfer silicone  
Standard packing 02525  
2 x 50 ml cartridges  
6 mixing cannulas, green

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