VOCO

IonoStar<sup>®</sup> Molar

(GB)

CE 0482

In accordance with DIN EN ISO 9917-1

Instructions for use

**Product description:**

IonoStar Molar is a radiopaque glass ionomer bulk filling material in user-friendly capsules for direct application. In addition to high compressive strength, it boasts good adhesion to enamel and dentine and releases fluoride ions.

IonoStar Molar is available in the shades A1, A2 and A3. The usable quantity of cement per capsule is at least 0.15 ml.

**Indications:**

- Restorations of non occlusion-bearing class I cavities
- Semi-permanent restorations of class I and II cavities
- Restorations of cervical lesions, class V cavities, root caries
- Restorations of class III cavities
- Restoration of deciduous teeth
- Base/liner
- Core build-up
- Temporary restorations

**Application:**

**Preparation:**

Prepare the cavity using a minimally invasive approach (removing only carious tooth substance). There is no need for undercuts. The layer thickness of the restorations must not be less than 1 mm. Avoid thinly tapering margins. Clean the cavity, rinse with water and dab or blow dry. Do not overdry the dentine: it should have a moist, shiny appearance. Place a matrix, if necessary. If the residual dentine thickness is insufficient, the area close to the pulp should be covered with a suitable calcium hydroxide product.

**Activating the application capsule and mixing the material:**

Only extract the application capsule from the tubular bag immediately before use. During activation, the liquid from the plunger must be pressed into the inside of the capsule as fully as possible. To achieve this, the end of application capsule where the plunger of liquid protrudes is slowly pressed down onto a hard surface (e.g. a table) as far as the plunger will move (see fig. 1). The capsule must only be activated directly before use and loaded immediately into a high-frequency mixer with approx. 4,000 oscillations per minute. Mix the activated application capsule for 10 - 15 s (mixing for longer reduces the working and setting times), remove it from the mixer without delay and insert it in the holder of an application device (e.g. AC Applicator, see corresponding instructions for use).

**Application:**

Carefully open the outlet by lifting the cannula until you feel a stop (see fig. 2), and apply the material directly into the cavity within 15 s of the end of mixing by compressing the lever on the applicator slowly and evenly. The cannula can be aligned as required by turning the whole application capsule inside the holder. Working time and consistency may be influenced by the respective type of capsule mixer and mixing time.

Fig. 1: Activating the capsule

Fig. 2: Opening the capsule

Following mixing, IonoStar Molar can be worked for at least 1.5 minutes (at room temperature). Prevent the material from being contaminated with moisture and from drying out. If it is not possible to keep the working area dry, a protective varnish should be applied after the contours/modelling have been finished (see corresponding instructions for use). The restoration can be finished, e.g. with a diamond bur or a polisher, 5 - 6 minutes after the start of mixing at the earliest. A protective varnish should then be applied again (see corresponding instructions for use).

Mixing time (23 °C)	10 - 15 s
End of mixing to application (23 °C)	15 s
Working time (23 °C)	1.5 min
Setting time (37 °C)	3 - 4 min
Total setting time	5 - 6 min

**Information, precautionary measures:**

IonoStar Molar contains polyacrylic acid. IonoStar Molar should not be used in case of known hypersensitivities (allergies) to this ingredient.

Do not apply directly to dentine near the pulp or to exposed pulp. In case of contact with the eyes, rinse immediately with plenty of water and consult a medical specialist. Rinse the instruments used immediately with cold water before residues of material cure. Ensure that the application capsule is mixed immediately after activation.

**Storage instructions and application method:**

Store IonoStar Molar in a dry place at temperatures of 15 °C - 23 °C. Only remove it from the foil immediately before use. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB2520 E1 0516 99 © by VOCO
**VOCO GmbH**