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**Bite registration material, silicone based, addition curing
fast setting – high final hardness**


Fig. 1

1. Mixing and Dosing

Extrude material from the cartridge with Automix2-system gun (Fig. 1). Place cartridge into mixing gun and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices (Fig. 2). Insert mixing cannula into guide grooves on cartridge and cannula and rotate clockwise to a stop (Fig. 3). Store cartridge with mixing cannula attached for subsequent use.


Fig. 2

2. Instruction for use

Apply milkbite directly from the mixing cannula, or using a contouring tip, onto the teeth (Fig. 4). Let the registration cure in the required bite position and remove it from the mouth after 50 seconds.

Fig. 3

3. Disinfection

After removal from the mouth rinse registrate under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde or other usual solutions.


Fig. 4

Important working hints

- Remove any retraction astringents or mouth wash solutions prior to the registration, these can influence the setting reaction.
- Do not work with condensation curing silicones.
- Cured silicone materials are chemically stable – avoid staining on clothing.
- Leave no residues of the material in the mouth.
- Please follow the instructions of the safety data sheet!


Other information:

Silicone impression materials are proven a million times over, unwanted effects cannot be expected in proper use. However, immune reactions, e.g. allergies or irritations cannot be excluded basically. In the case of doubt, we recommend performing an allergy test before application.

Fig. 5

Only for dental use by qualified personnel.

Indications

for use:

- Bite registration with crown and bridge impressions and for orthodontic indications, key for intraoral registration

Technical Data:

- DIN 13903
- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** base: cream catalyst: white
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 30 sec.*
- **Setting time:** approx. 1 min. 20 sec.*
- **Time in mouth:** approx. 50 sec.
- **Final hardness:** approx. 90 Shore A
- **Linear dimensional change:** approx. 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

milkbite
Standard packing **03380**
2 x 50 ml cartridges
8 mixing cannulas, pink
8 contouring tips

Eco packing **03381**
8 x 50 ml cartridges
32 mixing cannulas, pink

Mixing cannulas **02707**
pink, 48 pcs.

Contouring tips **02648**
12 pcs.

greenbite apple
Standard packing **02641**
2 x 50 ml cartridges
8 mixing cannulas, green
8 contouring tips

greenbite colour
Standard packing **02207**
2 x 50 ml cartridges
8 mixing cannulas, green
8 contouring tips

glassbite
Standard packing **02571**
1 x 50 ml cartridge
4 mixing cannulas, green
4 contouring tips

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