monoprint supra

DE

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monoprint supra

Instructions for use

Transfer silicone, addition curing, medium consistency - medium bodied – snap set effect


Fig. 1

Fig. 2

Fig. 3


Fig. 4


Fig. 5

1. Preparation of the impression

Depending on the impression technique, select a suitable impression tray (custom or stock metal or hard resin tray). With custom and hard resin trays, ensure that they are of adequate dimensions and maintain their shape. Using individual trays, mechanical retentions can be made in addition (e. g. perforation in the marginal area). Carry out preparations according to the implant system, position transfer caps and impression posts, and in addition e.g. open and fit the impression tray (Fig. 1). For optimal adhesion we recommend coating all trays with a thin layer of **trayloc® A** (Fig. 2). If there is no or too little retention to the transfer posts, adhesion can also be improved by applying a thin layer of **trayloc® A** to them. Let dry for about 3 – 5 min. Impression posts and transfer caps must be cleaned and dry. Check for possible residues from latex gloves on the surface.

2. Mixing and Dosing

The cartridge is extruded with the Automix2 system gun (Fig. 3). After inserting the cartridge into the mixing gun, remove the cartridge cap by twisting. Express a small amount of material until silicone is dispensed evenly from the two openings (Fig. 4). To attach the mixing cannula note the guides on the cannula and cartridge. Insert mixing tip and rotate clockwise to a stop (Fig. 5). Dispense the material with steady pressure. After use leave the mixing cannula on the cartridge until the next use.

3. Procedure by taking an impression

At the start put an adequate amount of **monoprint supra** in the impression tray. Then inject intraorally around the relevant areas. Intra-oral tips or application syringes can be used. Bring the filled tray into position in the patient's mouth with light pressure. Fix the impression without pressure until setting is complete.

4. Disinfection

After removing from the mouth, rinse the impression under running lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde or other usual solutions for 15 minutes.

5. Model casting

The impression should not be cast before 30 min. have elapsed. Dimensional stability is guaranteed for 14 days. Recommended model material includes class III or IV dental gypsums and commercially available model resins.

6. Cleaning of trays

After setting, material can be removed with a blunt instrument. Remaining film of **trayloc® A** is dissolved by immersion in acetone or isopropanol (be sure to use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.


Important working hints

- The application of **trayloc® A** improves distinctly the reliable fixation of transfer posts and retention aids.
- Solutions used before taking the impression (e.g. retraction solutions) can influence the setting reaction and must therefore be removed completely.
- Do not process with condensation curing silicones.
- Latex gloves and latex-contaminated surfaces can influence setting.
- Do not leave any impression material residues in the mouth.
- Cured impression materials are chemically resistant – avoid getting them on clothing.
- Please follow the instructions of the safety data sheet!

Further information:

Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

Only for dental use by qualified personnel.

Indications for use:

- Transfer impressions in implantology
 - a) Transfer impressions with closed trays
 - b) Transfer impressions with opened trays
 - c) Direct impressions
- Over- and fixation impressions in the combined technique
- Single-phase impressions

Technical Data:

DIN EN ISO 4823 - Type 2

- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 1:1
- **Colour code:** Base: yellow Catalyst: white
- **Final hardness:** approx. 70 Shore A
- **Mixing time:** omitted, (Automix2-system)
- **Working time:** approx. 2 min.*
- **Setting time:** approx. 4 min.*
- **Time in mouth:** approx. 2 min.
- **Strain in compression:** approx. 3.0 %
- **Recovery from deformation:** approx. 99.8 %
- **Linear dimensional change:** approx. 0.2 %
- **Application:** at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperature accelerate, decreased temperatures retard a. m. times.

Ordering information:

monoprint supra
Standard packing **02525**
2 x 50 ml cartridges
6 mixing cannulas, green

trayloc® A
17 ml bottle with brush **03398**
15 ml spray bottle **03392**

Mixing cannulas **02705**
green, 48 pcs.

Intra-Oral-Tips **02345**
white, 96 pcs.

Impression syringes **02030**
50 pcs.

DETAX

GmbH & Co. KG
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

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