# picobloc

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## Light-Curing Blockout Material

### Description

picobloc is a light-curing material for blocking out and preparing dental workpieces.

### Composition

Preparation from acrylic resins and initiators

### Areas of Application

- For sealing saw cuts
- Blocking out cervical margins
- Blocking out for thermo-forming of bleaching splints
- Covering during sandblasting
- Blocking out dies.

### Application

**Blocking out dies:** A clean and dry die is guaranteed for preparation. This may not be previously treated with varnishes or anything similar. For small undercuts or indentations, apply **picobloc gel** sparingly with the supplied dispensing tip. Light-cure the die for 2-5 minutes in halogen or UV devices, or for 90 seconds in stroboscopic devices.

**Varnishing dies:** After polymerisation, the prepared dies can continue to be treated as usual with die spacers and inter-space varnishes (self or light-curing).

**Laseal** by picodent, is especially suitable for die curing. Follow the directions for use by the particular manufacturer.

**Blocking out models for bleaching splints:** Fabricate a plaster model as usual from a high-strength dental stone (Implant-rock, picodent) and trim this in order to ensure optimal adaption of the splint. Use a cannula to apply a 0.5-1 mm thick layer of picobloc gel to the labial surfaces of the teeth to be bleached up to about 1 mm beneath the free gingiva. Cure for about 5 minutes in halogen or UV devices, or for 180 seconds in stroboscopic devices.

### Duplication/Removal of Inhibition Layer

If you need to duplicate, you can only do this with picobloc gel. Remove the inhibition layer with alcohol and let the model dry well. If die spacers were still used after blocking out with **picobloc**, please observe the specific manufacturer's instructions about this.

### Warning

**picobloc** blockout materials must exclusively be used outside the oral cavity. May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Gently wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

**Please note that the curing time is critically dependent on the emission performance of the lamps used. Correct maintenance of the devices and regular replacement of the light fixtures contribute to consistent curing. Also follow the particular device manufacturer's instructions.**

### Storage Instructions

Store at 4-25 °C. Keep away from direct sunlight. Always seal containers after use. Do not use the material after the expiry date.

### Delivery Units

**picobloc:** 3 x 3 g **picobloc gel** red, incl. 12 dispensing tips

**The product was developed for use in the dental field and must be used in accordance with the directions for use. The manufacturer accepts no liability for damage caused by the product being used in other ways. Furthermore, the user is obliged to test the material on his or her own responsibility for suitability and usability for the intended purposes before it is used, especially if this is not specified in the directions for use.**

### Note on Disposal

Dispose of contents/container in accordance with local regulation.

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picodent Dental-Produktions- und Vertriebs GmbH