DETAX

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Precision impression silicone for dental prosthetics

INDICATION

silaplast Futur, hydro C Putty, exaplast putty – tray material for correction impressions (two-step), double-mix technique (single-step) and anatomic impressions.

silasoft® Normal, silasoft® Special, hydro C Wash, exasoft wash – correction material, sprayable, for correction impressions (two-step), double-mix technique and anatomic impression (single-step), relining impression only for silasoft® Special and silasoft® Normal.

silone® – for taking impressions of dentate jaws, anatomic impressions, onlay and inlay impressions.

sta-seal – for taking impressions of edentulous jaws, mucostatic functional and relining impressions.

sta-seal f – for functional impression-taking with extended curing phase, functional margin design.

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

Do not use for the treatment of pregnant/breastfeeding women and children (catalyst component contains alkyl stannosiloxanes).

INTENDED USERS

Dentist

PRODUCT DESCRIPTION

The CDS silicones are condensation-curing, elastomeric impression materials.

The CDS silicones have differing consistencies in accordance with DIN EN ISO 4823 type 0 – 3 (see table on page 16).

METERING AND MIXING (see table on page 16)

silaplast Futur – with catp universal FUTUR, paste / with silasoft® catf, liquid

silone®, sta-seal – with catp universal FUTUR, paste / with silasoft® catf, liquid

sta-seal f – with sta-seal catf, liquid

hydro C Putty, hydro C Wash – with hydro C Activator, paste

exaplast putty, exasoft wash – with exagel cat / with silasoft® catf, liquid

silaplast Futur, hydro C Putty, exaplast putty, sta-seal f – manual metering and kneading

silasoft® Normal, silasoft® Special, silone®, sta-seal, hydro C Wash, exasoft wash – manual metering and mixing with a spatula

PROCEDURE

TRAY, ADHESIVE

Select a suitable impression tray in accordance with the impression technique. To optimize silicone adhesion, apply sili Adhesive Fluid or spray it on (Fig. 1).

MANUAL METERING AND MIXING

Putty – The metering is carried out with the enclosed metering spoon (Fig. 2a) or with the dosing syringe (Fig. 2b).

Paste – The metering is carried out from the tube onto a mixing block in the required strand length (Fig. 3).

Catalyst liquid – The metering is carried out with the enclosed metering spoon or by dripping, respectively.

Mixing is done manually by kneading or using a spatula on a mixing block until a homogeneous mixture or uniform color, respectively, is achieved. The mixing time is 30 seconds (Fig. 4). **silasoft® N / silasoft® S / hydro C Wash / exasoft wash** can be taken up and applied using a disposable or an impression syringe (Fig. 5a/5b). For **sta-seal f**, use of an easy Application Syringe is recommended for individual margin design (Fig. 5c). After metering, immediately close the container again.

DISINFECTION

After taking the impression out of the mouth, rinse it under running lukewarm water. Subsequent disinfection can be carried out with suitable commercially available disinfecting solutions in accordance with the manufacturer's instructions.

PRODUCING THE LABORATORY MODEL

With dental plasters type III – V, this can be done already about 30 minutes, up to no later than 72 h, after removal from the mouth.

TRAY CLEANING

Remove cured material mechanically. Remove sili Adhesive Fluid with sili Adhesive Solvent. Use only in well-ventilated areas. Then clean and disinfect the tray as usual.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment during the finishing of the cured material.
- ▶ Avoid eye contact with catalyst, rinse eyes immediately with plenty of flowing water if necessary, consult ophthalmologist immediately.
- ▶ Cured impression materials are chemically resistant – avoid stains on clothing.
- ▶ Do not leave any residues of the impression material in the patient's mouth or sulcus, and make sure to avoid inadvertent pressing of material into areas outside the oral cavity through mouth-antrum connections or malformations of the palate.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ The good flowability of **silasoft® N** guarantees excellent detail fidelity. Because of this very light-bodied consistency, in rare cases **silasoft® N** is prone to partial segregation. Before metering the catalyst, homogenize any segregated **silasoft® N** with a spatula on the mixing block.
- ▶ **exasoft wash / hydroC Wash**: Any silicone oil segregated by storage can be mixed normally with the other components and has no effect on processing. The curing behavior is not affected thereby.
- ▶ Overmetering of catalyst and higher temperatures accelerate, undermetering and lower temperatures delay the curing.
- ▶ Carefully close the catalyst bottle or tube immediately after use.
- ▶ Read and understand the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority of the Member State in which the user and/or patient is established.

CONTRAINDICATION






Do not use this material in case of known allergies to any of the ingredients, or contact allergies. Silicone impression materials have been tried and tested millions of times; undesired effects are not to be expected in case of proper application. However, as a matter of principle immune reactions such as allergies or irritations cannot be ruled out. In case of doubt, we recommend performing an allergy test before application.

SIDE EFFECTS


Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Technical Specifications

	DIN EN ISO 4823 TYPE	Dosage Base mass	Dosage cat f	Dosage cat p	Product color Base/Catalyst	Open time* min. approx.	Time in mouth min. approx.	Curing time* min. approx.
silaplast Futur	type 0, putty	1 ×  = 13.5 g	1 ×  = 0.5 g		white / light blue	1:15	2:30	3:45
				6 cm  = 0.3 g		white / blue	1:15	3:30
silasoft® Normal	type 3, light-bodied	10 cm  = 5 g	10 ×  = 0.2 g		pink opaque / light blue	1:15	2:30	3:45
				10 cm  = 0.5 g		pink opaque / blue	1:15	2:30
silasoft® Special	type 3, light-bodied	10 cm  = 4.5 g	8 ×  = 0.16 g		pink-transp. / light blue	1:30	2:30	4:00
				10 cm  = 0.5 g		pink-transp. / blue	1:30	2:30
silone®	type 2, medium-bodied	1 cm  = 1 g	1 ×  = 0.02 g		russet / blue	1:15	2:30	4:00
				1 cm  = 0.05 g		russet / blue	1:15	2:30
sta-seal	type 2, medium-bodied	1 cm  = 1.1 g	1 ×  = 0.02 g		lime green / blue	1:15	2:30	3:45
				1 cm  = 0.05 g		lime green / blue	1:15	2:30
sta-seal f	type 1, heavy-bodied	1 cm  = 1.7 g	1 ×  = 0.02 g		green / blue	2:00	4:00	6:00
hydro C Putty	type 0, putty	1 ×  = 13.5 g		8 cm  = 0.4 g	pastel green / blue	1:15	3:30	4:45
hydro C Wash	type 3, light-bodied	10 cm  = 5 g		10 cm  = 0.5 g	pink opaque / blue	1:15	2:45	4:00
exaplast putty	type 0, putty	1 ×  = 13.5 g		6 cm  = 0.3 g	yellow / blue	1:15	3:30	4:45
exasoft wash	type 3, light-bodied	10 cm  = 5 g		10 cm  = 0.5 g	pink / blue	1:15	2:30	3:45

Processing: At 23 °C ± 2 °C,
50 ± 5% relative humidity

Storage:

* from start of mixing at 23 °C ± 2 °C,
50 ± 5 % rel. humidity. Higher temperatures shorten, lower temperatures extend the indicated times.

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Ordering Information

silaplast Futur <u>Standard packing</u> 1 jar 900 ml	02001	silasoft® Special <u>Standard packing</u> 1 tube 160 ml	02275	hydro C <u>Set</u> Putty, jar 1.000 ml Wash, tube 150 ml Activator, tube 60 ml	02869	<u>Standard packing cat</u> 1 tube 60 ml	03125
<u>4-pack</u> 4 jars each 900 ml	02002	<u>4-pack</u> 4 tubes each 160 ml	02276	<u>Standard packing Putty</u> 1 jar 1.000 ml	02242	Catalysts <u>cat p universal Futur, paste</u> 5 x 35 ml suitable for silaplast, silasoft® Normal, silasoft® Special, silone®, sta-seal	02015
<u>Clinic packing</u> 1 jar 5.400 ml	02003	silone® <u>Standard packing</u> 1 tube 160 ml	02008	<u>Eco packing Putty</u> 1 jar 5.400 ml	02471	<u>silaplast Futur cat f, liquid</u> 5 x 50 ml suitable for silaplast	02016
<u>Refill bag</u> 1 jar 900 ml	02681	<u>3-pack</u> 3 tubes each 160 ml	02009	<u>Standard packing Wash</u> 1 tube 150 ml	02243	<u>silasoft® cat f, liquid</u> 5 x 10 ml suitable for silasoft® Normal, silasoft® Special	02017
silasoft® Normal <u>Standard packing</u> 1 tube 160 ml	02004	sta-seal <u>Standard packing</u> 1 tube 160 ml	02010	<u>Standard packing Activator</u> 1 tube 60 ml	02245	<u>sta-seal cat f, liquid</u> 5 x 10 ml suitable for silone®, sta-seal, sta-seal f	02018
<u>4-pack</u> 4 tubes each 160 ml	02005	<u>3-pack</u> 3 tubes each 160 ml	02011	exaplast <u>Set</u> putty, jar 910 ml wash, tube 140 ml cat, tube 60 ml	03133		
<u>Standard packing bottles</u> 1 bottle 160 ml	02320	sta-seal f <u>Standard packing</u> 3 measuring syringes each 80 ml 1 sta-seal cat f bottle, 10 ml	02012	<u>Standard packing putty</u> 1 jar 910 ml	02907		
<u>4-pack</u> 4 bottles each 160 ml	02321	<u>Large packing</u> 9 measuring syringes each 80 ml 3 sta-seal cat f bottles, 10 ml	02013	<u>Standard packing wash</u> 1 tube 140 ml	02910		

DETAX

GmbH & Co. KG

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

Telefon: 07243/510-0 · Fax: 07243/510-100

www.detax.de · post@detax.de

Made in
Germany

11/2020
3.1

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).