# silasoft® direct

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**GB**

Instructions for use

**Precision impression material, C-silicone based,  
light bodied consistency – normal setting, 10:1**

Fig. 1

## 1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **sili Adhesive Spray** or **sili Adhesive liquid** onto the tray and let dry for 60 sec (Fig. 1). Take preliminary impression with **silaplast FUTUR**. To facilitate the flowing out of excessive light body material from the impression cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended.

Fig. 2

## 2. Mixing and Dosing

Extrusion of the impression material from the cartridges by means of Automix2 gun (Fig. 2). Place cartridge into mixing gun, remove cartridge plug by rotation. Eject a small amount of material (bleeding) to control that equal amounts are being extruded from both orifices (Fig. 3). Insert mixing cannula into guide grooves on the cartridge and rotate anti-clockwise to a stop. Squeeze mixing gun handle for individual dosage (Fig. 4). For application into sulcus area use Intra-Oral tips yellow – coloured or a disposable syringe (Fig. 5).


Fig. 3

## 3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.


Fig. 4

## 4. Model casting

Not before 30 min. after removal from the mouth, at the latest after 24 hours. Recommended materials: Dental plasters of class III and IV-also commercially available model acrylics.


Fig. 5

## 5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

## 6. Cleaning of trays

Remove set material with a blunt instrument. **sili Spray** or **sili liquid** film can be easily removed with **sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean and disinfect as usual.

## Important working hints


- Increased temperatures accelerate, decreased temperatures retard the setting time.
- After use reclose cartridge orifices carefully.
- Do not intermix with addition curing silicones.
- Avoid staining of clothing. Cured impression materials cannot be removed.
- Remove all residues of impression material from the patient's mouth.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek medical care.
- Please follow the instructions of the safety data sheet!

Avoid contact with eyes and skin!  
Wear suitable eye protection and gloves!

Other information:

Silicone impression materials are proven a million times over, unwanted effects cannot be expected in proper use. However, immune reactions, e.g. allergies or irritations cannot be excluded basically. In the case of doubt, we recommend performing an allergy test before application.

Only for dental use by qualified personnel.

## Indications for use:

- Preliminary impressions for the double impression or double mix technique (wash)
- Reline impressions

## Technical Data:

DIN EN ISO 4823 – Type 3

- **Mixed volume:** 50 ml (cartridge)
- **Mixing ratio:** 10:1
- **Colour code:** base: violet catalyst: blue
- **Mixing time:** omitted (Automix2-system)
- **Working time:** approx. 1 min. 15 sec.\*
- **Setting time:** approx. 4 min. 15 sec.\*
- **Time in mouth:** approx. 3 min.
- **Strain in compression:** approx. 5 %
- **Recovery from deformation:** approx. 98.5 %
- **Linear dimensional change:** approx. 0.5 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

\* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

## Ordering information:

- silasoft® direct Automix2**
- |                                 |              |
|---------------------------------|--------------|
| <b>Standard packing</b>         | <b>02809</b> |
| 2 x 50 ml cartridges            |              |
| 12 mixing cannulas, yellow 10:1 |              |
| <b>Eco packing</b>              | <b>02810</b> |
| 8 x 50 ml cartridges            |              |
| 48 mixing cannulas, yellow 10:1 |              |
| 48 Intra-Oral-Tips, yellow      |              |
| <b>Mixing cannulas</b>          | <b>02807</b> |
| 48 pcs., yellow 10:1            |              |
| <b>Intra Oral Tips</b>          | <b>02712</b> |
| 96 pcs., yellow                 |              |
| <b>Mixing gun Automix2</b>      | <b>02808</b> |
| 1 pc. 10:1                      |              |

**DETAX**

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 0 72 43/510-0 · Fax: 0 72 43/510-100  
www.detax.de · post@detax.de

Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

Made in Germany

07/2017

# "C" silicone system

putty,  
heavy bodied

kneadable,  
very high  
consistency

## silaplast FUTUR

high final hardness  
normal setting

Preliminary impression  
material for corrections of  
preparations, tray material,  
overimpressions, wax-bite

Wash  
material  
light bodied

low viscosity  
low con-  
sistency

## silasoft® N

normal setting

Correction impression  
material for the double  
impression and double  
mix technique, relinings

Wash  
material  
light bodied

low viscosity  
low con-  
sistency

## silasoft® S

normal setting

Correction impression  
material for the double  
impression and double  
mix technique, relinings

Functional  
heavy bodied

high con-  
sistency and  
plasticity

## sta-seal f

medium final hardness  
plastic phase  
extended setting time

Functional impressions,  
functional border moulding

Catalyst  
paste

Catalyst  
liquid

## cat<sup>p</sup> / cat<sup>f</sup> FUTUR

inert system  
non-irritating

cat<sup>p</sup> universal paste,  
cat<sup>f</sup> silaplast liquid,  
silasoft® cat<sup>f</sup>, liquid,

Special  
instrument

Accessory

## Deta-Cut

## sili adhesive

## sili solvent

Sterilizable, curved, for  
cutting the preliminary  
impression

adhesive for C-silicones  
removal of the adhesive  
and cleaning of trays

All DETAX condensation curing silicones are portioned with dispensing aids (tray, measuring bowl, measuring scale etc.) and kneaded by hand or mixed homogeneously on a mixing pad.

### Ordering information:

**silaplast FUTUR**  
Standard packing **02001**  
jar, 900 ml

4-pack **02002**  
4 x 900 ml

clinic packing **02003**  
1 x 5400 ml

refill bag **02681**  
900 ml

**silasoft® N**  
standard  
tube-packing **02004**  
160 ml

4-pack tubes **02005**  
4 x 160 ml

standard  
bottle-packing **02320**  
160 ml

4-pack bottles **02321**  
4 x 160 ml

**silasoft® S**  
standard  
tube-packing **02275**  
160 ml

4-pack tubes **02276**  
4 x 160 ml

**sta-seal f**  
standard-packing **02012**  
measuring syringes  
3 x 80 ml

catalyst liquid bottle 10 ml  
large-packing **02013**  
measuring syringes 9 x 80 ml  
catalyst liquid bottles  
3 x 10 ml

**cat<sup>p</sup>/cat<sup>f</sup>**  
cat<sup>p</sup> universal FUTUR  
paste, 5 x 35 ml **02015**

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