Deutsch
English
Español
Français
Italiano

Monopolarer Elektrodenhandgriff smartPen

Monopolar electrode handle smartPen

Mango portaelectrodos monopolaes smartPen

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Contents

1	Notices Concerning this Document	34
1.1	Safety.....	34
1.2	Terms & Acronyms.....	34
1.3	Validity of this Document.....	35
1.4	Symbols Used in this Document	36
2	Product Liability and Warranty.....	36
2.1	General Information	36
2.2	Hotline.....	37
2.3	User’s Inspection.....	37
2.4	Warranty	37
3	Intended Use	38
3.1	Intended Purpose.....	38
3.2	Ambient conditions.....	38
3.3	Users.....	38
3.4	Approved Rated Voltages	38
3.5	Warnings	40
4	Use	43
4.1	Before First Use and Before Each Further Use	43
4.2	Application	43
4.2.1	Inserting the Active Electrode.....	43
4.2.2	Connecting to the Electrosurgery Unit	45
4.2.3	Activate HF current	45
4.3	After Completion of the Application.....	47
5	Processing	48
5.1	Basic Information on Processing.....	48
5.2	Reprocessing Limitation and Restrictions.....	49
5.3	Preparation for Cleaning.....	49
5.4	Manual prewash	50
5.5	Machine Cleaning and Disinfection	51
5.5.1	Phases of the machine cleaning process	52
5.5.2	Drying.....	52
5.6	Sterilization	53
5.6.1	Before Sterilization.....	53
5.6.2	Additional Guidelines.....	54

5.7	Storage and Transportation	54
5.8	Functional Check	55
6	Disposal	56
6.1	Packaging	56
6.2	Disposal of the Product	56
6.3	National Regulations	56

1 Notices Concerning this Document

1.1 Safety

Failure to observe this document can lead to serious or even lethal injury to the patient or the user!

Improper handling and care as well as non-intended use can lead to premature wear and/or pose a risk to patients and users!

It is the operator's responsibility to ensure that all personnel handling the product have understood and do observe the notes and instructions in this document.

- Every user is required to read this document completely and follow it carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- The Instructions for Use provided in this document are not intended as a replacement for the Operating Instructions electro surgery unit, the foot switch and the neutral electrode. In particular the safety notes therein must always be observed!

1.2 Terms & Acronyms

Designation	Description
CJD	Creutzfeldt-Jakob disease
HF	High Frequency
NE	Neutral electrode
NMR	Nuclear Magnetic Resonance
NTP	Low temperature plasma sterilization
MRT	Magnetic Resonance Tomography
RDG	Cleaning and disinfection device

1.3 Validity of this Document

This document applies to the following products:

Designation	Connection	REF
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	KLS Martin	80-217-13-04
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	3-Pin	80-210-13-04
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	Erbe VIO, ICC, ACC	80-210-18-04
smartPen, electrode handle, 1 key, Ø 4 mm, cable length 4 m	KLS Martin	80-217-05-04
smartPen, electrode handle, without button, Ø 4 mm, cable length 4 m	KLS Martin	80-220-02-04
smartPen, electrode handle, without button Ø 4 mm, cable length 4 m	Erbe VIO, ICC, ACC	80-220-31-04
smartPen, electrode handle, without button, Ø 4 mm, cable length 4 m	Bovie	80-220-41-04
smartPen, electrode handle, 2 keys, Ø2,4 mm, cable length 4 m	3-Pin	80-210-21-04
smartPen, electrode handle, rocker switch, Ø 2,4 mm, cable length 4 m	3-Pin	80-210-23-04
smartPen, electrode handle, 2 keys, Ø 4 mm	without cable	80-227-00-04
Monopolar cable for smartPen, REF 80-227-00-04	KLS Martin	80-227-01-04
Monopolar cable for smartPen, REF 80-227-00-04	Erbe VIO, ICC, ACC	80-227-02-04
Monopolar cable for smartPen, REF 80-227-00-04	3-pin	80-227-03-04
smartPen, electrode handle, 1 key, Ø 1,6 mm, cable length 4 m	KLS Martin	80-214-13-04
smartPen, electrode handle, without button, Ø 1,6 mm, cable length 4 m	KLS Martin	80-214-12-04

1.4 Symbols Used in this Document

Words marked with ® are trademarks owned by the respective proprietors.

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

Life hazard or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

Risk of injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device/machine failure, etc.)!

2 Product Liability and Warranty

2.1 General Information

We thank you for having decided to buy one of our products. This product bears the CE-marking, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

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info@klsmartin.com · www.klsmartin.com

2.2 Hotline

- Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact the Product Management:

Tel: +49 7461 706-243

Fax: +49 7461 706-190

- Should you have any technical questions or questions concerning maintenance contracts or training courses, please contact our Martin Service Center:

Tel: +49 7461 706-343

Fax: +49 7461 706-484

E-mail: service@klsmartin.com

NOTICE

Each product and its packaging are usually marked with a specific lot number (LOT) and a reference number (REF). Please always indicate LOT and REF in case of a complaint.

2.3 User's Inspection

- Immediately upon receipt, the goods must be checked for completeness and potential damage in transit.
- Notice of any such damage must be given immediately.
- The smartPen is packed in a fully hygienic but non-sterile packaging.
- In case of return shipments, send only clean and disinfected products in sterile packing.

2.4 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form.

3 Intended Use

3.1 Intended Purpose

Handheld monopolar electrode handle designed for the connection to an electrosurgical device and enabling monopolar cutting and coagulation on a patient by the use of an active electrode.

3.2 Ambient conditions

Store the smartPen at room temperature in a sterilization container/soft package in a clean and dry environment with controlled humidity. Information regarding the minimum and maximum ambient temperature for storage and transport, see section 5.7 "Storage and Transportation", page 54.

- Protect against direct sunlight. Storage in the original packaging up to the first use is recommended.

3.3 Users

The smartPen may be used only by competent and trained medical staff who are familiar with the safety precautions and risks of electrosurgery.

The smartPen may be used only by persons specially instructed and trained in its use.

3.4 Approved Rated Voltages

Risk of injury from electric shock!

There is risk of electric shock if the output voltage of the electrosurgery unit is higher than the rated voltage of the smartPen.

- Make sure that the output voltage of the electrosurgery unit cannot exceed the rated voltage of the smartPen!

The power ranges of active electrodes depend on the surgical application. Selection of the correct performance parameters is the user's responsibility.

Bezeichnung	Connection	Rated voltage [kVp]	REF
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	KLS Martin	6,3	80-217-13-04
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	3-Pin	6,3	80-210-13-04
smartPen, electrode handle, 2 keys, Ø 4 mm, cable length 4 m	Erbe VIO, ICC, ACC	6,3	80-210-18-04
smartPen, electrode handle, 1 key, Ø 4 mm, cable length 4 m	KLS Martin	6,3	80-217-05-04
smartPen, electrode handle, without button, Ø 4 mm, cable length 4 m	KLS Martin	6,3	80-220-02-04
smartPen, electrode handle, without button, Ø 4 mm, cable length 4 m	Erbe VIO, ICC, ACC	6,3	80-220-31-04
smartPen, electrode handle, without button, Ø 4 mm, cable length 4 m	Bovie	6,3	80-220-41-04
smartPen, electrode handle, 2 keys, Ø 2,4 mm, cable length 4 m	3-Pin	6,3	80-210-21-04
smartPen, electrode handle, rocker switch, Ø 2,4 mm, cable length 4 m	3-Pin	6,3	80-210-23-04
smartPen, electrode handle, 2 keys, Ø 4 mm	ohne Kabel	6,3	80-227-00-04
Monopolar cable for smartPen REF 80-227-00-04	KLS Martin	6,3	80-227-01-04
Monopolar cable for smartPen, REF 80-227-00-04	Erbe VIO, ICC, ACC	6,3	80-227-02-04
Monopolar cable for smartPen, REF 80-227-00-04	3-pin	6,3	80-227-03-04
smartPen, electrode handle, 1 key, Ø 1,6 mm, cable length 4 m	KLS Martin	6,3	80-214-13-04
smartPen, electrode handle, without button, Ø 1,6 mm, cable length 4 m	KLS Martin	6,3	80-214-12-04

3.5 Warnings

 WARNING**Risk of injury due to electric shock or severe burns!**

- **Active implants such as pacemakers or implanted electrodes**
can be damaged or disturbed by application of electrosurgery. Moreover, implants can change the current density, leading to unnoticed endogenous burns.
- Consult a cardiologist before the operation!
- **Metallic shunts in the HF current path**
- The smartPen electrode handle may be used only with HF electrodes by Gebrüder Martin.
 - Compatible electrodes are listed in the HF accessories catalog by KLS Martin
 - Compatible electrodes can be found on the KLS Martin website

The patient should not come into contact with metallic parts that are grounded or have significant electrical capacitances to the ground.

- We recommend using antistatic cloths for this purpose. In cases of doubt, feel free to contact the Martin Service Center for further information, see section 2.2 “Hotline”, page 37.
 - Skin-to-skin contact, e.g. between the arms and the body of the patient, must be avoided, for example by placing dry gauze in between.
 - **Alternating HF current paths**
In case of direct skin contact of cables, capacitive currents may cause burns.
 - Lay out the supply cable to the smartPen in the operating area so that it comes into contact neither with the patient nor with other cables.
 - **Uncontrolled emission of HF current!**
 - Avoid accidental activation of HF current by unintentional activation of a key, rocker switch or foot switch!
 - Keep magnetic materials and magnetic fields (especially MRI, NMR spectroscopy) away from the foot switch!
-

 WARNING**Peril to life or danger of serious physical injury from prion infection!**

Destroy products immediately after any potential contact with prions and do not use again under any circumstances. There is risk of contamination with Creutzfeldt-Jakob Disease (CJD)!

 **WARNING**

Risk of burns from depositing surgical instruments on the patient!

Uncontrolled emission of HF current may occur if the smartPen is placed on the patient and thus comes into contact with tissue.

- Do not deposit a temporarily unused smartPen on the patient! Avoid accidental contact with tissue as well as contact with low-resistance materials, for example liquids!
 - Coagulate or cut only if the contact surface of the active electrode of the smartPen is visible. Avoid touching other metallic instruments or objects!
-

 **CAUTION**

Risk of injury from tissue damage!

HF energy can interfere with physiological monitoring devices used together with electrosurgery on the patient.

- Attach monitoring electrodes as far away as possible from HF electrodes.
- Do not use needle electrodes for monitoring.

During procedures where HF current flows through any body part with small cross section, **tissue damage** can occur! For example, this applies to bony parts with minimally conductive upper skin layer. To avoid unintentional coagulation, the application of bipolar technology may be an option.

- Do not apply over bony edges and scar tissue.
 - For current flow, select parts between the active electrode and neutral electrode with good conductivity.
-

 **WARNING**

Risk of fire or explosion from flammable or explosive substances!

Materials saturated with oxygen, such as cotton, gauze or insulating blankets, may catch fire from the sparks that occur during the intended use of the electrosurgical device.

If e.g. during monopolar coagulation, the active electrode is only slightly in contact with the tissue, by design sparking between the tissue and the electrode will be observed. Sparks in the vicinity of flammable gases (e.g. explosive anesthetic gases) can lead to explosions.

There is the additional risk of accumulation of flammable liquids under the patient or in anatomical cavities and recesses.

- Never use the instruments in the presence of combustible or explosive substances!
-

⚠ CAUTION**Risk of burns or malfunction from improper application!**

The function and product service life of the smartPen can be affected by improper handling. Cables with damaged or brittle insulation or cable with broken electrical lines may cause fire and burns to the user or patient.

- Protect the smartPen from mechanical damage, do not throw or drop, do not apply force!
- Connect the smartPen only to electrosurgery units intended for this type.
- The smartPen, made of composite materials, cannot be repaired. In case of impaired function, please contact the Martin Service Center, see section 2.2 "Hotline", page 37.
- Regularly check the cables for damage, when required also using a magnifying glass.

Improper handling of the cable may damage the extremely thin and fine wires and the passivation layers that protect them, and can thus lead to premature malfunction.

Therefore, avoid:

- Pulling on the cable.
 - Squeezing and bending the lines/cables.
 - Tight wrapping of the cables.
 - Wrapping the cable around the smartPen.
 - Handling of cables entails danger of tripping!
 - Store cables in such a way that the insulation of the cable is protected against mechanical damage.
-

4 Use

4.1 Before First Use and Before Each Further Use

WARNING

Risk of infection in case of non-sterile handling!

The smartPen is supplied in **non-sterile** packaging.

- Clean, disinfect and sterilize before first use and before each further use, see section 5 “Process”, page 48.
-

- Please note the safety instructions, see section 3.5 “Warnings”, page 40.
- Inspect the smartPen for damage and completeness. Especially check the insulation of the cable for damage before each use.
- Check “Approved Rated Voltages”, of the active electrode, see section 3.4, page 38!
- Perform functional test, see section 5.8 “Functional Check”, page 55.

4.2 Application

4.2.1 Inserting the Active Electrode

WARNING

Risk of injury from electric shock and/or severe burns!


- The smartPen electrode handle may be used only with HF electrodes by Gebrüder Martin.
 - Disconnect smartPen from the electrosurgery unit before inserting the active electrode!
 - The used active electrode can be very hot after the application!
 - Do not touch hot surfaces!
-

NOTICE


Possible damage to the active electrode!

Active electrodes must not be exposed to leverage and are to be protected against mechanical stress!

- Insert a sterile active electrode suitable for the surgical application (with a connection diameter of 4.0 mm, 2.4 mm, or 1.6 mm, depending on the electrode handle model) into the electrode socket of the reprocessed smartPen until the end stop is reached. The dimensions and corresponding tolerances for the electrode socket intended for connection of active electrodes to the smartPen can be provided upon request.
 - The active electrode can be inserted secured against rotation in various positions (in case of electrodes with hexagon).
 - Insert the electrode into the electrode carrier of the smartPen in a straight direction; jamming (as well as rotating movement) can lead to damage due to the leverage effect.
 - The hexagon of the active electrode or the end of the insulation, respectively, must be fully covered by the tip of the smartPen.
 - No metal parts may be exposed or visible, respectively, at the transition area from electrode to handle. The insulation of an insulated active electrode must be inserted into the electrode carrier down to the end stop.
 - The following illustrations show a properly inserted electrode.

- A handle-electrode combination where metal parts (see marking in the following illustration) are exposed or visible, respectively, **must not be used**, as in the worst case this can result in a shunt and thus cause burns.

In case of questions, please contact Gebrüder Martin.

4.2.2 Connecting to the Electrosurgery Unit

The connection to the Gebrüder Martin electrosurgery unit is established via the cable on the smartPen or a separate connection cable.

Risk of injury from electric shock and/or severe burns!

- Connect the connection cable with smartPen only to the monopolar socket on a deactivated electrosurgery unit or an electrosurgery unit in **STANDBY** mode.
 - Connect the smartPen only to an electrosurgery unit suitable for that type.
-
- Before inserting the connection cable, ensure that the connector matches the monopolar socket on the electrosurgery unit.
 - If a connector requires a specific orientation relative to the socket before insertion, ensure that the orientation is correct. Forceful insertion results in damage to the connector and/or the socket!
 - Fully insert the connector into the electrosurgery unit until the mechanical limit stop is reached. The plug connection must provide for proper, permanent contact.
 - Turn on electrosurgery unit and select cutting or coagulation current so that the rated voltage of the active electrode is not exceeded, see section 3.4 “Approved Rated Voltages”, page 38.

4.2.3 Activate HF current

The HF current can be activated using the keys on the smartPen. Generally, the smartPen can be activated using a foot switch as well.

smartPen with 2 keys

- The yellow key is used for activation of the current type Monopolar Cutting.
- The blue key is used for activation of the current type Monopolar Coagulation.

smartPen with 1 key

- The key is used for activation of the current type pre-selected on the electrosurgery unit.

smartPen without keys

- The respective HF current is activated using the foot switch.

smartPen with rocker switch

- The yellow rocker switch is used for activation of the current type Monopolar Cutting.
- The blue rocker switch is used for activation of the current type Monopolar Coagulation.

⚠ WARNING**Risk of burns from high HF current concentration!**

Unplugging the connector of the connection cable of the smartPen and plugging it into another device socket during application can result in unintentional energy emission due to different settings of the two device sockets and thus cause burns.

- Observe the specific default settings of the respective device sockets of the electrosurgery unit!
-

⚠ WARNING**Danger of burns from power modifications!**

If the coagulation or cutting effect is reduced during operation, never increase the output power of the electrosurgery unit without prior check!

- Specifically, check the following:
 - Cleanliness and wear of the distal end of the active electrode.
 - Proper placement of the neutral electrode over the entire contact surface on the patient's skin.
 - Flawless contact of all HF plugs and connection cables.
 - Function of the activation switches on the smartPen.
 - Function of the foot switch.
 - Foot switch or insulation of the HF components (external damage).
-

⚠ WARNING**Risk of burns from depositing surgical instruments on the patient!**

Uncontrolled emission of HF current may occur if the smartPen is placed on the patient and thus comes into contact with tissue.

- Do not deposit a temporarily unused smartPen on the patient! Avoid accidental contact with tissue as well as contact with low-resistance materials, for example liquids!
 - Coagulate or cut only if the contact surface of the active electrode of the smartPen is visible. Avoid touching other metallic instruments or objects!
-

4.3 After Completion of the Application

Risk of burns due to hot surface!

The active electrode of the smartPen can be very hot after the application!

- Do not touch hot surfaces!
-
- Deactivate electrosurgery unit.
 - Unplug the connector from the monopolar socket of the electrosurgery unit.
 - Do not pull on the cable when loosening the plug connection. Pulling on the cable may damage the extremely thin and fine wires and the passivation layers that protect them, and can thus lead to premature malfunction.
 - Disconnect the active electrode from the smartPen.

5 Processing

5.1 Basic Information on Processing

The smartPen is supplied in non-sterile condition.

- Sterilize the smartPen before initial Use, see section 5.6 “Sterilization”, page 53.
-

NOTICE

Risk of damage!

To prevent tissue residues from drying, we recommend reprocessing of the contaminated instrument **immediately** after use.

The operator/user is responsible for cleaning, disinfecting and sterilizing the instrument. It is essential that national regulations and restrictions be observed.

Use appropriate and approved sterile packing (e.g. in accordance with EN 868, ISO 11607) for sterilization, subsequent transport and storage.

The parameters for cleaning and disinfecting the product are to be determined by the operator/user, e.g. as part of the hygiene plan. The treatment agents must be selected to be compatible with the components, as well as hygienically and microbiologically effective.

A thorough final rinse must be performed to ensure that no residues of cleaning or care agents are left on any of the individual components after the reprocessing.

NOTICE

- Observe all applicable national regulations in regard to reprocessing of the products in case of patients suffering from Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of the disease.
 - Note that successful reprocessing of this medical product can be guaranteed only after prior validation of the reprocessing process. The responsibility for this lies with the operator/reprocessor.
 - In the interest of process capability and better and more reliable cleaning results, machine reprocessing should always be preferred to manual cleaning, as mechanical cleaning processes can be validated.
 - Due to process tolerances, the data provided by the manufacturer serve only as reference values for the evaluation of the reprocessing processes carried out by the operator/reprocessor.
-

5.2 Reprocessing Limitation and Restrictions

With each application, a progression of the aging of the product is associated. Product life is usually determined by wear and damage during use. This wear is technical and inevitable. Therefore, the product must undergo a thorough visual and functional check before each application.

For information on functional testing, see chapter 5.8 Functional Testing.

Non-functional products must be replaced.

- Dispose of the smartPen according to the national regulations, see section 6 “Disposal”, page 56.

Peril to life or danger of serious physical injury from prion infection!

- Observe the applicable national regulations regarding the reprocessing of the products in case of patients with Creutzfeldt-Jakob Disease (CJD), suspected CJD or any variants thereof.

5.3 Preparation for Cleaning

Danger of serious injury from electric shock!

Prior to cleaning or disinfection, be sure to de-energize the system completely:

- Disconnect smartPen from electrosurgery unit.
- Immediately after use in the operating theater, wipe the smartPen clean using a synthetic fleece soaked with surface disinfecting agent.

Risk of burns due to hot surface!

Active electrodes can be very hot after application!

- Do not touch hot surfaces!
-
- Pull active electrode out of electrode socket.

5.4 Manual prewash

NOTICE

Risk of damage from tissue residues!

Dried or fixated surgical residues may damage the surfaces of the smartPen.

- The smartPen must be thoroughly cleaned and disinfected **immediately** after each use.
 - Perform the cleaning process as early as possible (ideally within 2 h).
 - Do not use fixating pre-cleaning or fixating disinfecting agents.

Attached particles or secretions may render later cleaning or sterilization difficult or impossible and can therefore result in malfunction.

Step 1

- Clean the handle surface and the contaminated portion of the cable with a soft brush (e.g. insitumed REF Med. 100.33) under cold running water until no more dirt is visible.
- Disconnect the active electrode from the handle. When a smartPen REF 80-227-00-04 is used, unplug cable.

Step 2

- Immerse the smartPen into a cold water bath (<35°C/<95°F) for 5 min. All surfaces should be covered completely.
- During this time, the activation buttons (if present, depending on the model) should be pressed three times.
- Clean the surfaces again with a soft brush (e.g. insitumed REF Med. 100.33).
- Flush the lumen of the electrode connection with (e.g. with a syringe 10 ml and a cannula held at an angle of 45° to the lumen).

Step 3

- Ultrasound cleaning with e.g. Cidezyme (ASP, Johnson & Johnson), working solution 0.8% at max. 40°C (104°F) for 5 min. Observe the specific information by the manufacturer.

Step 4

- Rinse the electrode of the smartPen (e.g. with a water jet gun for 10 s, distance of 6–8 cm).
- Rinse the handle for 1 min under cold running tap water.
- Check all parts for visible contamination.
- If there is any residual dirt, repeat the process.

NOTICE

Risk of damage in case of unsuitable cleaning agents!

- Unsuitable cleaning agents or excessive temperatures may result in damage to the instrument.
 - Use only cleaning and disinfecting agents that have been approved for use with plastics and stainless steel.
 - Use only approved cleaning and disinfecting agents according to the manufacturer's instructions.
 - Observe the information regarding concentration, temperature and exposure time.
 - Use of hydrogen peroxide can lead to premature wear of the instrument and thus reduce its lifetime.
-

5.5 Machine Cleaning and Disinfection

Recommendation for machine cleaning: Validation in accordance with EN ISO 15883 and accompanying information in accordance with EN ISO 17664. Cleaning with different process parameters, validated acc. to EN ISO 15883, is possible.

Like other surgical instruments, the smartPen can be reprocessed by machine. Manual pre-cleaning is required only in cases of particularly heavy contamination, see section 5.4 "Manual prewash", page 50.

NOTICE

- The tray baskets selected for the cleaning process must be loaded with such a quantity and in such a manner that no impairment of the cleaning results is to be expected.
 - Do not clean the smartPen together with sharp-edged or pointed objects.
 - Even after only a short period, surfaces with metallic coating can show material changes upon use of an alkaline cleaning agent and/or an acidic neutralizer. This may lead to impairment of functionality. The material is best protected by use of a pH-neutral or mildly alkaline cleaning agent/treatment agent.
-

5.5.1 Phases of the machine cleaning process

- Place the smartPen, including the cable, into the RDG load carrier (e.g. Miele disinfecter G 7836 CD).

Phase	Step	Temperature	Time	Water quality	Chemistry
1	Pre-cleaning	Cold water max. 40°C (104°F)	2 min	Drinking water	--
2	Cleaning	55°C (131°F), max. 65°C (149°F)	5 min	Fully demineralized water if possible	pH-neutral to alkaline, e.g. neodisher® Medizym Working solution 0.8%
3	Neutraliza- tion	Cold water max. 40°C (104°F)	3 min	Fully demineralized water if possible	--
4	Interim rinsing	Cold water max. 40°C (104°F)	2 min	Fully demineralized water if possible	--
5	Thermal disinfection	e.g. 90°C (194°F) A0 3000, see EN ISO 15883	5 min	Fully demineralized water if possible	--
6	Drying	--	According to program	--	--

Instrument care agents or surface tension-reducing agents can be added if they are approved for the medical instruments and are a part of the validated program.

- After this procedure, all parts should be checked for visible dirt. If there is any residual dirt, repeat the process.

CAUTION

Danger of burns caused by disinfection and cleaning agents!

Insufficiently removed disinfection agents can lead to burns.

- Securely remove any residues of disinfection agents by through rinsing.

5.5.2 Drying

NOTICE

Ensure reliable drying!

Reliable drying is an essential factor for successful sterilization!

5.6 Sterilization

5.6.1 Before Sterilization

NOTICE

- Ensure that the smartPen has been cleaned.
- Intactness of the insulation, cleanliness and integrity of the instrument must be visually ascertained.
- The active electrode must have been removed from the smartPen!
- Repeat the cleaning if a visual check reveals tissue residues on the smartPen, or if the active electrode is still in the smartPen, see section 5.4 "Manual prewash", page 50.
- Properly pack the smartPen according to ISO 11607 and EN 868 in order to maintain the sterile barrier after removal from the sterilizer.
- Please follow the recommendations of the manufacturer of the sterilizer as regards handling, loading and parameter setting.
- Sterilization at high temperatures or long sterilization times reduces the lifespan of the instruments.

⚠ WARNING

Risk of infection in case of non-sterile handling!

Improper sterilization and non-sterile handling of the instruments can entail severe health risks for the patient.

- The active electrodes must be properly protected so that the sterile barrier is not damaged.
- In case of return shipments, send only clean and disinfected products in sterile packing.
- Always sterilize the active electrode and smartPen separately.
- Sterilization must be done using a validated steam sterilization method, e.g. in a sterilizer according to EN 285:2009 and validated according to ISO 176651:2006.
- The operator must ensure that sterility is maintained after the sterilization process.
- Use appropriate and approved sterilization packaging (e.g. in accordance with EN 868, ISO 11607) or sterilization containers for storage and transport.

⚠ WARNING

Risk of infection and risk of damage in case of impermissible sterilization!

Impermissible sterilization can lead to severe health risks for the patient and damage to the smartPen.

- Do not sterilize the smartPen with hot air or flash autoclaving methods.
- Do not use any plasma sterilization methods (NTP Sterilization)!

! WARNING

Peril to life or danger of serious physical injury from prion infection!

- Destroy the smartPen immediately after potential contact with prions (risk of CJD contamination!) and do not use it again under any circumstances!

The temperature in the autoclave must not exceed 138°C (280.4°F) and 18 min holding time, since otherwise the insulation and other plastic parts might be damaged.

Steam sterilization with fractionated pre-vacuum	
Minimum requirement	3 min at 132°C (269.6°F)
Maximum sterilization temperature	138°C (280.4°F)
Maximum sterilization time	18 min
Drying time	according to validation/program

Other times and temperatures can also be used as long as they are within the appropriate range.

The operator must ensure that sterility is maintained after the sterilization process.

Follow the handling and loading recommendations by the manufacturer of the sterilization device!

5.6.2 Additional Guidelines

Cleaning, care and sterilization instructions and other regulations are available in the Red Brochure entitled “Proper Maintenance of Instruments”, published by the AKI (“Arbeitskreis Instrumentenaufbereitung” – Instrument Processing Working Group), which can be obtained from Gebrüder Martin GmbH & Co. KG upon request. Useful information on instrument processing is also available via our website at: <http://kls-martin.com/phnet/>.

5.7 Storage and Transportation

Information on minimum and maximum temperature for storage and transportation:

Designation	Monopolar electrode handle smartPen	
Environmental conditions for storage and for transport	Ambient temperature	-20°C to +50°C (-4°F à +122°F)
	Relative humidity (no condensing moisture!)	10–75%
	Air pressure	500–1,600 mbar

- Store the smartPen in a clean, cool and dry place.
- Protect smartPen against mechanical damage.
- When the variant smartPen REF 80-227-00-04 is used, unplug cable.
 - Store cables in such a way that the insulation is protected against mechanical damage.
- Store and transport the smartPen in safe containers or sterile packing. Use appropriate certified sterilization containers, e.g. according to EN 868, ISO 11607.

5.8 Functional Check

Risk of injury from instruments or components that have not been checked!

Non-functional or defective instruments and components can endanger the patient or user, as well as affect the intended function of the instrument.

- Always keep electrosurgical instruments and their components in flawless working order.
 - Never use electrosurgical instruments and their components if they are damaged by improper handling or transport. If instruments or components are damaged, they must be disposed of immediately.
 - Use of components or accessories by other manufacturers may constitute a source of danger.
-

Before each application, perform the following functional checks:

- Check the “Approved Rated Voltages ” of the active electrode, see section 3.4 “Approved Rated Voltages”, page 38!
- Visual check for potential damage (check for cracks, fractures or deformation):
 - Electrode handle
 - Electrode socket
 - Active electrode
 - Insulation of the connection cable
 - Plug connections
- Check for easy movement of all control elements, especially keys or rocker switches and foot switches.
- Check the function of the activation switch (e.g. by activation).

6 Disposal


6.1 Packaging

Gebrüder Martin will, as a matter of course, take back the full packaging if so desired, to recycle as many parts of the packaging as possible.

If you do not wish to make use of this offer, you can dispose of the packaging with the normal paper or domestic garbage.

6.2 Disposal of the Product

In designing the product, we tried to avoid using composite materials wherever possible. This allows a high degree of recycling. We therefore offer to take the product back for proper disposal and recycling.

	<p>Marking of electric and electronic equipment in accordance with Directive 2002/96/EC (WEEE Directive) and the German Electrical and Electronic Equipment Act (ElektroG)</p> <p>This symbol on the product or its packaging indicates that the product may not be disposed of as normal household garbage.</p>
---	---

6.3 National Regulations

The national regulations and disposal provisions must be observed for all disposal measures.

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Martin Nederland/Marned B.V.
KLS Martin UK Ltd.