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Long-term temporary crown & bridge material, self curing, bis-acrylic-composite 10:1


Fig. 1

1. Preparation of the impression

Before preparing the stump or bridge preparation or a planned extraction, take an anatomic impression using alginate or silicone. Interdental ridges should be excised to increase the amount of material available for finishing the temporary (Fig. 1). If there are gaps in the molar region, a connecting bar between prepared teeth is created by cutting a groove between the abutment teeth in the impression. In case of gaps in the anterior region prior to impression-taking, use prefabricated teeth as placeholders. Block these solidly, e.g. with wax. Alginate impressions should be kept at 100% rel. humidity until filling with **tempofit® semi**.


Fig. 2

2. Use of the dispensing gun

Cartridge material is extruded using the Automix2-System dispensing gun (Fig. 2). After inserting cartridge into dispensing gun, unscrew cap and carefully extrude a small amount of material to insure proper flow from both ducts (Fig. 3). Attach mixing cannula by aligning notched side to cartridge and turning clockwise until fastened (Fig. 4). Dispenser now ready to mix silicone in any amount needed. Leave mixing cannula on cartridge after use. Use only blue **tempofit® semi** 10:1 blue mixing cannulas!


Fig. 3

3. Mixing and application

tempofit® semi base material and catalyst are automatically mixed in the correct proportion when extruding by dispensing gun. Before beginning work with the Automix2 dispensing gun, always extrude a small (pea-sized) amount and discard (Fig. 5). Apply **tempofit® semi** to deepest parts of impression first and evenly fill impression up to gingival level (Fig. 6).


Fig. 4

4. Intraoral modelling of the temporary

Seat filled impression onto prepared teeth and remove any excess material with a plastic tool. **tempofit® semi** will attain a hard-elastic consistency after 2 - 3 min. after extrusion, at which time temporary and impression can be removed from mouth. Be sure to closely follow the progress of polymerisation, as removal is only possible while material is still elastic (Fig. 7).


Fig. 5

5. Hardening and finishing

Remove temporary from impression and remove excess material and proximal undercuts. After 6 - 7 min., material is completely cured and hard enough to be trimmed and polished with finishing burs. Contact with oxygen creates an inhibition layer on the temporary, which should be removed using ethyl or isopropyl alcohol. Instead of polishing, the finished temporary can be surface-sealed with **tempofit® finish**.


Fig. 6

6. Fastening the temporary

The temporaries can be fixed using commercially available temporary dental cements (eugenol-free) and, in case of longer wear times, using permanent dental cements. Remove all excess cement after incorporation. Fully eliminate cement residues when removing the temporary.


Fig. 7


7. Individualization and repair

Perform shade individualization, reinforcement or repair using **smartrepair®**. Esthetic adaptation of the gingival areas is possible using **smartrepair® rosa**. Thoroughly rinse and dry the temporary. Roughen the area to be treated and be sure to remove any residual dust from grinding. Apply **smartbond®** outside the mouth to cover all areas and rub in using **smartbrush**. Disperse with compressed air and allow to dry for 1 min. Apply **smartrepair®** and light-cure from all sides for at least 40 sec. Finish the light-curing process outside the mouth by light-curing from all sides for additional 30 sec. Trim the repaired areas and polish the temporary or seal with **tempofit® finish**.

Important hints

- **tempofit® semi** can be shade matched using the **smartrepair®** system and composites.
- **tempofit® semi** can be used to fabricate long-term temporary restorations in the anterior and posterior regions and for single crowns and bridges up to four units with one pontic. The minimum occlusal thickness is 1.5 mm, circular 0.8 mm. The connector cross-section for pontics in the anterior region is at least 12 mm², in the posterior region at least 15 mm². Failure to observe these specifications can lead to an undesirable outcome.
- To ensure a homogeneous mixing: Prior to insert the mixing cannula, **always** eject a small amount of material to ensure extruding **equal amounts** of base and catalyst **from both orifices** (Fig. 3).
- Avoid contact with skin and eyes. In the event of accidental contact, rinse with plenty of running water, and consult an specialist if necessary.
- Please follow the instructions of the safety data sheet!

For use by trained specialists.

Indications for use:


- Long term temporaries
- Temporary crowns and bridges

Contra-indication:

- **tempofit® semi** based on (meth)acrylates. In the case of hypersensitivities (allergies) to these ingredients, sensitization may be caused. Do not use **tempofit® semi** with persons having such a predisposition.

Technical data:

- **Mixed volume:** 50 ml (cartridges)
- **Mixing ratio:** 10:1
- **Product colours:** A2 / A3 / A3.5 (as indicated on package)
- **Dosage and application time:** approx. 45 sec.*
- **Setting time in the mouth:** approx. 2 - 3 min.* (hard-elastic phase, from beginning of mixing)
- **Setting time:** approx. 6 - 7 min.* (until setting is completed)
- **Application:** At 23 °C ± 2 °C / 73 °F ± 3.6 °F
- **Storage:**

* At 23 °C: Increased temperatures accelerate, decreased temperatures prolong these times.

Physical properties:

- **Compressive strength:** > 250 MPa
- **Flexural strength:** > 80 MPa
- **Diametral tensile strength:** > 40 MPa
- **Barcol hardness:** > 37

Ordering information:

tempofit® semi 10:1	
Standard pack	
Cartridge of 50 ml, 10:1	
T-Mixer 10:1, blue	
A2	04052
A3	04053
A3.5	04054

Dispensing gun	
10:1	02808

Mixing cannulas	
10:1	04055
T-Mixer 10:1, blue 25 pcs.	

tempofit® finish	
Presentación	
normal	03536
botella	4 ml
Pincel desechable	10 pzs.

smartrepair®	
smartrepair® set	
application syringe	02694
A2	1,5 g
application syringe	
A3,5	1,5 g
smartbond® bottle	5 ml
smartbrush	5 pcs.
application needles	6 pcs.
flowchart	1 pc.

Standard packing	
application syringe	1,5 g
application needles	5 pcs.
A2	02772
A3	02862
A3.5	02695

smartrepair® rosa	
application syringe	03364
application needles	3 g
	8 pcs.

smartbond®	
bottle	02696
smartbrush	5 ml
	5 pcs.

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