

**BioSonic<sup>®</sup>**

Ultrasonic Cleaning Systems

[biosonic.coltene.com](http://biosonic.coltene.com)

 **COLTENE**

2018



BioSonic UC150



BioSonic UC125



BioSonic UC50DB



BioSonic UC95



BioSonic UC300



BioSonic UC300R



BioSonic UC100



BioSonic UC1

1988

# BioSonic®

## An essential part of your Infection Control System

Ultrasonic cleaning is a process created by high frequency soundwaves. The soundwaves create high energy cavitation. During cavitation millions of tiny bubbles form and then implode, releasing enormous amounts of energy and shock waves. This powerful scouring action reaches into minute crevices which cannot be reached by manual scrubbing. The combination of energy and specially formulated solutions make ultrasonic cleaning the most effective method for removing gross and microscopic debris.

BioSonic is a complete ultrasonic cleaning system which consists of three different cleaner unit sizes, a full line of cleaning solutions and unit accessories which provide you with the flexibility to efficiently clean and disinfect instruments, prostheses and other items.

### BioSonic units feature unique product design and optimize time and safety in the office or laboratory

#### Gentle but effective cleaning

- Instruments get cleaned thoroughly but stay in shape and sharp even after frequent cleaning cycles

#### Quiet performance

- Closed lid operation
- Plastic housing

#### User safety

- Built in drain with extra long drain tube for safe disposal of contaminated solutions (UC125 and UC300)
- Raised ridges inside baskets for safe removal of instruments

#### Long lasting attractive appearance

- Wear-resistant plastic housing which will not stain or dent

\*Some units shown have been discontinued



# BioSonic® UC50DB



## Technical Information

Tank capacity	3.2 liters
Fill capacity	2.0 liters
Overall Dimensions (l x w x h)	27.3 x 21.9 x 20.3 cm
Tank Dimensions (l x w x d)	19.7 x 14.6 x 10.5 cm
Output	50 watts
Operating Frequencies	53 kHz

## It's a small giant!

### Small on the outside, big on the inside

- The compact BioSonic UC50DB fits any work surface but its 2 litres capacity accommodates double ended instruments or two 600 ml beakers. The FingerGuard™ plastic basket with metal bottom provides superior wave transmission and prevents dulling of instruments.

### Convenient Chairside Placement

- The unique design of cover and housing minimizes noise and prevents the escape of airborne contaminants enabling chairside cleaning of patients' dental appliances using the beaker option.



## BioSonic® ST16 Dental Soaking Tray



### Technical Information

Tank capacity	0.47 liters
Dimensions	6.4 x 22.9 x 12.7 cm

The BioSonic Instrument Soaking Tray is a sure way to reduce the risks associated with instrument handling prior to sterilization. Pre-soaking dental instruments prior to ultrasonic cleaning is the clinically approved method for effective infection control, while reducing tedious hand scrubbing.

- Compact design is ideal for operatory
- Chemically resistant inner tray lifts out automatically when the lid is opened
- Contemporary design matches many office decors



# BioSonic® UC125




## Technical Information

Tank capacity	4.8 liters
Fill capacity	3.8 liters
Overall Dimensions (l x w x h)	42 x 26.7 x 35 cm
Tank Dimensions (l x w x d)	28.6 x 14.3 x 12.7 cm
Output	100 watts
Operating Frequencies	55 kHz *

\* C-Power – power is automatically adjusted in the tank based on load requirements.


## Customized operation and reliability at its best




**Multiple Languages**  
Choose one of six different languages



**Degas Button**  
Easily prepare your solution for maximum effective cleaning




**Selectable Run Times**  
Customize cleaning procedures of your office by selecting run times for different tank loads



**Countdown Timer**  
Know exactly when your instruments will be ready for the autoclave




**Solution Tracking**  
See exactly how long your solution has been in use




**Solution heating**  
Preheat your solution to desired temperature

## Cross Flow Baskets

- Enhanced Instrument Cleaning
- Easy Lift-out Handles



**UC151XD**  
Full Size Basket



**UC152XD**  
Three Quarter Size Basket

# BioSonic® UC300



## Technical Information

Tank capacity	14.4 liters
Fill capacity	11 liters
Overall Dimensions (l × w × h)	48.2 × 30.5 × 31.1 cm
Tank Dimensions (l × w × d)	35.6 × 22.9 × 17.1 cm
Output	300 watts
Operating Frequencies	40 kHz

## A large capacity cleaner for your practice



### Cycle large volumes of instruments

- Twin-baskets provide uninterrupted instrument processing
- Clean two baskets of instruments as you process third basket for sterilization and fourth basket for pre-soak

### Greater cleaning efficiency

- Large basket allows you to clean 3 to 4 times as many instruments
- Basket dividers facilitate sorting and minimize handling of instruments
- Easily holds larger items (e.g. film processor rollers and articulators)

### Clean instruments and prostheses – simultaneously

- Clean both your instruments and dental appliances – simultaneously!



# BioSonic® Solutions

The solution for your cleaning needs

## MeterDose® Bottles – accurate solution dosage

Dispense the exact amount needed. Simply squeeze the desired amount into the MeterDose® measuring chamber; tilt and pour.

- No Waste
- No Spills
- No Mess
- Easy To Use
- Economical



## BioSonic® UC30, 473 ml

Super Concentrate General Purpose

**Concentration** 1.6% to be used in the ultrasonic bath diluted yield 30 liters.

**Indication:** UC30 is a non-ionic detergent and is used for general cleaning. It removes most organic soils like blood, saliva and debris.

**Items:** Instruments crowns, bridges, dentures and ultrasonic scaler tips

**Special notes:** Use as a carrier bath in ultrasonic tank. Can also be used as instrument pre-soak and carrier bath in tank when using beakers. After cleaning thoroughly clean with tap water.



## BioSonic® UC31, Single Dose 20 × 15 ml

General Purpose + Super Rust Inhibitor

**Concentration** 0.4% to be used in the ultrasonic bath diluted yield 76 liters.

**Indication:** UC31 combines the proven cleaning action of a BioSonic solution with corrosion protection in the ultrasonic bath. It's used for general cleaning.

**Items:** Instruments crowns, bridges, burs, dentures and ultrasonic scaler tips

**Special notes:** Use as a carrier bath in ultrasonic tank. Can also be used as instrument pre-soak and carrier bath in tank when using beakers. After cleaning thoroughly clean with tap water.





## BioSonic® UC32, 236 ml Enzymatic Cleaner

**Concentration** 0.8% to be used in the ultrasonic bath diluted yield 30 liters.

**Indication:** UC32 combines the non-ionic detergency of the UC30 with enzymes that act to break down stubborn dirt particles. Effectively breaks down and loosens debris hardened onto instruments. Minimizes hand scrubbing.

**Items:** Instruments, crowns, bridges, dentures and ultrasonic scaler tips

**Special notes:** Use as a carrier bath in ultrasonic tank. Can also be used as instrument pre-soak and carrier bath in tank when using beakers. After cleaning thoroughly clean with tap water.



## BioSonic® UC34, 840 g Plaster and Stone Remover

**Concentration** 1 scoop (tbs) per 100 ml of solution to be used in the ultrasonic bath diluted yield 11.4 liters.

**Indication:** Concentrated powder for easy removal of plaster and stone. Built in scoop for measuring correct dosage.

**Items:** Attachments, articulators, instruments, castings and items that can not be sandblasted

**Special notes:** Do not use with aluminium. Can also be used directly in ultrasonic tank. After cleaning thoroughly rinse with tap water. Replace after each use.



## BioSonic® UC38, 236 ml Germicidal Cleaner

**Concentration** 0.8% to be used in the ultrasonic bath diluted yield 30 liters.

**Indication:** Removes blood, debris, saliva and kills HIV virus and tuberculosis while it cleans instruments in the ultrasonic bath.

**Items:** Instruments, crowns, bridges and ultrasonic scaler tips

**Special notes:** Do not clean dentures or plastics. Use as a carrier bath in ultrasonic tank. Can also be used as instrument pre-soak and carrier bath in tank when using beakers. After cleaning thoroughly clean with tap water.



## BioSonic® UC39, 473 ml Cement Remover

**Concentration** 40.0% to be used in the ultrasonic bath diluted yield 32 uses.

**Indication:** Cleaner with cement softener. Removal of temporary and permanent cements (including zinc-phosphate, polycarboxylate).

**Items:** Crowns, veneers, brackets and instruments

**Special notes:** Do not use with aluminium. Do not use directly in ultrasonic tank. After cleaning thoroughly rinse with tap water.

# BioSonic® Ultrasonic Cleaning Systems



## Order Information

### UC300

<b>UC300</b>	<b>BioSonic® Ultrasonic Cleaner</b>	
	Tank Cover	1 pc.
	Enzymatic Cleaner (UC32)	1 pc.
	Owner's guide	1 pc.

#### Accessories for UC300

<b>UC310</b>	Large FingerGuard™ Instrument Basket (incl. 5 basket section inserts)	1 pc.
<b>UC311</b>	Twin FingerGuard™ Instrument Basket (incl. 2 basket section inserts)	1 pc.
<b>UC330</b>	2-Beaker Positioning Cover	1 pc.
<b>UC335</b>	Basket Section Inserts	2 pcs.
<b>UC340</b>	Replacement Tank Cover	1 pc.
<b>UC350</b>	Replacement Drain Filter	3 pcs.
<b>UC355</b>	Replacement Drain Tube	1 pc.
<b>UC325</b>	Cassette Basket	1 pc.

#### Accessories for all units

<b>ST16</b>	Soaking Tray	1 pc.
<b>UCB2</b>	Basket for small parts	1 pc.
<b>UC53</b>	600 ml Beaker, Cover and Positioning Ring	1 pc.

### UC125

<b>UC125H</b>	<b>BioSonic Ultrasonic Cleaner with heating</b>	
	Tank Cover	1 pc.
	General Purpose + Super Rust Inhibitor (UC31)	20 × 15 ml
	Owner's guide	1 pc.

#### UC125 BioSonic® Ultrasonic Cleaner without heating

	Tank Cover	1 pc.
	General Purpose + Super Rust Inhibitor (UC31)	20 × 15 ml
	Owner's guide	1 pc.

#### Accessories for UC125

<b>UC155XD</b>	Accessory Kit (¾ FingerGuard™ Instrument Basket for UC125, 2 Beakers, Covers, Positioning Rings and 2 One-Beaker Positioning Racks)	
<b>UC151XD</b>	Full Size FingerGuard™ Instrument Basket for UC125	1 pc.
<b>UC152XD</b>	Three-Quarter FingerGuard™ Instrument Basket for UC125	1 pc.
<b>UC153</b>	1 Beaker Positioning Rack	1 pc.
<b>UC350</b>	Replacement Drain Filter	3 pcs.
<b>UC355</b>	Replacement Drain Tube	1 pc.

### UC50DB

<b>UC50DB</b>	<b>BioSonic Ultrasonic Cleaner</b>	
	Tank Cover	1 pc.
	FingerGuard™ Instrument Basket	1 pc.
	Owner's guide	1 pc.


#### Accessories for UC50DB

<b>UC59</b>	Accessory Kit (FingerGuard™ Instrument Basket, 600 ml Beaker, Cover and Positioning Ring, 1 Beaker Positioning Rack, Basket for small parts)	
<b>UC58</b>	Plastic FingerGuard™ Instrument Basket	
<b>UC153</b>	1 Beaker Positioning Rack	1 pc.

#### BioSonic Cleaning Solutions

<b>UC30</b>	Super Concentrate General Purpose Solution in MeterDose® Bottle	473 ml
<b>UC31</b>	General Purpose + Super Rust Inhibitor Single Dose	20 × 15 ml
<b>UC32</b>	Enzymatic Cleaner MeterDose® Bottle	236 ml
<b>UC34</b>	Plaster & Stone Remover Powder Concentrate in Jar	840 g
<b>UC38</b>	Germicidal Cleaner MeterDose® Bottle	236 ml
<b>UC39</b>	Cement Remover MeterDose® Bottle	473 ml

© COLTENE – www.coltene.com



#### Coltène/Whaledent Ltd.

The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF/UK  
T +44 1444 235486  
F +44 1444 870640  
info.uk@coltene.com

#### Coltène/Whaledent AG

Feldwiesenstrasse 20  
9450 Altstätten/Switzerland  
T +41 71 757 53 00  
F +41 71 757 53 01  
info.ch@coltene.com

#### Coltène/Whaledent Inc.

235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223/USA  
T +1 330 916 8800  
F +1 330 916 7077  
info.us@coltene.com

COLTENE