








# FenderMate® Prime

Protection and Matrix  
for primary teeth.

**DIRECTA**  
directadental.com





Use rubber dam or secure with waxed dental floss.

### FenderMate® Prime

**long**

t= 0,10 mm



**short**

t= 0,08 mm



602815 FenderMate Prime Assortment 2x18 pcs  
602816 FenderMate Prime Short, neon yellow 18 pcs  
602817 FenderMate Prime Long, neon green 18 pcs  
602808 FenderMate Prime Value Pack, Short, neon yellow 100 pcs  
602809 FenderMate Prime Value Pack, Long, neon green 100 pcs

**REF**

**Rx Only**



FenderMate® Prime (Registered Design and Patent pending) is a registered trademark owned by Directa.



**DIRECTA AB**

P.O. Box 723, Finvids väg 8  
SE-194 27 Upplands Väsby,  
Sweden

Tel: +46 8 506 505 75

**DIRECTA INC**

64 Barnabas Road, Unit 3  
Newtown, CT 06470,  
United States

203.491.2273

106876-2019-04-12