



Expert information

REMOT

CE

Manufacturer

lege artis Pharma GmbH + Co. KG
P. O. Box 60, D-72132 Dettenhausen
Breitwasenring 1, D-72135 Dettenhausen
Telephone +49 (0) 71 57 / 56 45 - 0
Fax +49 (0) 71 57 / 56 45 50
E-Mail: info@legeartis.de
www.legeartis.de

Name

REMOT prophylaxis paste
For removal of dental plaque, prophylaxis of caries and final polish.
For dental use.

Composition

Zirconium(IV)-oxide, silica, cocamidopropyl betaine, macrogols, sorbitol, glycerol, chlorhexidine digluconate and peppermint oil

Indications

REMOT removes plaque and therefore prevents caries and periodontal disease.
Also suitable for final polish.

Notice

The well-balanced combination of abrasive materials effects a careful removal of the dental plaque without damaging the enamel.

Treatment is not recommended in the case of known allergies to one of the ingredients.

Mode of application

REMOT will be applied with rubber cup and prophy brushes.

Shelf life

REMOT shall not be used after the expiry date.

Administrative form and package size

75 g prophylaxis paste Item number 0032317

Date of revision

2019-04

Symbol of "Expiry date"



Symbol of "Batch number"



Symbol of "Manufacturer"



Symbol of "Catalogue number"

Shows Item number of the manufacturer

