

Hands-Free Gingival Retraction System**SECTION I - Product/ Company Identification**

Product Trade Name: GingiTrac™ Retraction

Part (Item) Number: REF 360100, REF 360101, REF 360102

Manufacturer: CENTRIX INC

Address: 770 RIVER ROAD, SHELTON, CT 06484

Emergency Telephone Number: Toll-Free (800)235-5862, (203)929-5582

SECTION II - Hazardous Ingredients/Identity Information

Chemical Family: Condensation Silicone

Hazardous Components: Contains silicic acid esters and dibutyltin dilaurate

SECTION III - Hazards Identification

HAZARD CLASSIFICATION: Irritant

SECTION IV - First Aid Measures

Skin: Wash with plenty of water and soap.

Ingestion: Consult a physician.

Inhalation: If affected, remove to fresh air.

Eyes: Flush immediately with plenty of water. If irritation persists, see eye doctor.

SECTION V - Fire-Fighting Measures

Suitable extinguishing media: Waterspray, foam, powder, carbon dioxide

Unusual fire and explosion hazards: Material can burn generating irritating gases

Specific intervention methods: Appropriate protective equipment must be worn in case of fire

SECTION VI - Accidental Release Measures

Recovery: Contain and absorb spilled material using an appropriate spill absorbent. Collect in a suitable container for disposal. Cleaning/Decontamination: Wipe up any affected area with absorbent material.

Disposal: Dispose of material in accordance with local, state, and federal requirements.

SECTION VII - Handling and Storage

Handling: Does not require specific technical measures. Note that handling with latex or rubber gloves may inhibit set.

Use vinyl gloves. Air sensitive material; keep container tightly closed when not in use.

Storage Containers: Store in supplied packaging (keep container tightly closed)

Storage Incompatibilities: Strong oxidizing agents

Storage Other precautions: Store away from heat sources. Temperature will affect set time.

SECTION VIII - Exposure Controls / Personal Protection

Eye protection: Eye protection recommended

Hands protection: Impervious gloves are recommended

Other equipment: None required

SECTION IX - Physical and Chemical Properties

Appearance: Green base/white catalyst

Vapor pressure:< 5 mm Hg @ 25 C

Odor: Slight sweet odor

Water solubility: Insoluble in water

Boiling point: >150°C Specific gravity:> 1

Flash point: >300°F (material will ash)

Volatile part by weight: 30 percent

SECTION X - Stability and Reactivity

Stability: Stable

Materials to avoid: No dangerous reactions known

Hazardous decomposition products: Silicic acid esters may hydrolyze to form 2-methoxyethanol which can pose additional health hazards

SECTION XI - Toxicological Information

Acute toxicity: Non-toxic

Eye contact: Eye irritant

SECTION XII - Ecological Information

Environmental fate: No data available

Degradability: Non-biodegradable

Ecotoxicity: Impact on aquatic environment: No data available. No affect anticipated.

SECTION XIII - Disposal Considerations

Dispose in accordance with federal, state and local regulations

SECTION XIV - Transport Information

DOT- not regulated

SECTION XV - Regulatory Information

Federal - SARA title III section 313: This product does not contain any substance subject to reporting requirements at or above deminimis quantities.

Federal - TSCA: All components are either listed or otherwise comply with requirements

International - DSL (Canada): All ingredients are listed on the Canadian Domestic Substance List (DSL)

International - WHMIS (Canada): Not a "controlled product" under the Canadian Workplace Hazardous Materials Information System- (WHMIS)

SECTION XVI - Other Information

None.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.



Making Dentistry Easier.SM

CENTRIX INCORPORATED

770 River Road, Shelton, CT 06484 USA
(800) 235-5862 Tel: (203) 929-5582
Fax: (203) 929-6804



CENTRIX EUROPE

Hauptstrasse 64 B, 65719 Hofheim, Germany