

Maximum power.
Whisper quiet.
Optimal comfort.

KaVo MASTERtorque LUX

- Precise performance with maximum power of up to 23 watts (M9000 L)/19 watts (M8700 L)
- Optimal comfort with Silence Technology
- Better safety & hygiene with Direct Stop Technology (DST)



KaVo MASTERtorque
LUX M9000 L

KaVo MASTERtorque
Mini LUX M8700 L



KaVo MASTER Series Turbines

KaVo MASTERtorque LUX M9000 L The Intelligent Turbine.

Power:

- 23 watts power
- 20% more torque and force
- Stable and efficient cutting

Direct Stop Technology (DST):

- Active, intelligent break system
- Bur stops in one second

Zero retraction:

- No suck-back of contaminated material
- Clean as a whistle – Longer lifetime

Silence Technology:

- 57 dB(A)
- Intelligent sound design eliminates high frequencies
- Pleasant run-down noise



KaVo MASTERtorque Mini LUX M8700 L Small Head. Safe Access.

Great visibility:

- 12.1 mm head size
- KaVo patented 100°/19° head angle
- 22° sight angle
- Compatible with long, standard and short burs

Power:

- 19 watts power
- Stable and efficient cutting

Direct Stop Technology (DST):

- Active, intelligent break system
- Bur stops in one second

Zero retraction:

- No suck-back of contaminated material
- Clean as a whistle – Longer lifetime

Silence Technology:

- 59 dB(A)
- Intelligent sound design eliminates high frequencies
- Pleasant run-down noise




NEW!

Intelligent features for your comfort:

- Plasmatec surface – Better grip
- Water filter easily exchangeable by dentist – Consistent spray performance
- Passive pressure regulation, TU-pressure from 2.1 to 3.5 bar – Consistent power

Proven KaVo quality features:

- Chuck retention up to 32 N – Safety and long life
 - Tungsten carbide guiding bush – Centric preparation and longer lifetime
 - Original KaVo ceramic bearings
 - 24 months warranty
-  Made in Germany

KaVo MASTERtorque Mini Kids Edition



MASTERtorque Mini LUX M8700 LK with funny dragon design "Koko"

Includes:

- 5 Booklets "Koko the humming dragon"
- 300 Sticker "Dragon Power"
- 1 Poster



KaVo master series

Article	Description	Material number
M9000 L	MASTERtorque LUX	1.008.7900
M9000 LS	MASTERtorque LUX	1.008.5400
M9000 L COLOR	MASTERtorque LUX (chocolate brown)	1.007.9900
M9000 L COLOR	MASTERtorque LUX (anthracite)	1.007.7100
M8700 L	MASTERtorque Mini LUX, for MULTIflex coupling	3.001.0000
M8700 LS	MASTERtorque Mini LUX, for Sirona coupling	3.001.0600
M8700 LK	MASTERtorque Mini Kids Edition	1.011.2330
MULTIflex 465 LED	MULTIflex LUX Turbine Attachment 465 LED	1.007.3197
QUATTROcare PLUS 2124 A	Handpiece Maintenance System	1.008.3805



KaVo. Dental Excellence.