
















#whdentalwerk
f @ in ▶
video.wh.com

Proxeo Aura
Air Polishing
System



Application	Setting	Powder	Spray head	
 <p>Supragingival</p>	 <p>Prophy</p>		 <p>Standard</p>	 <p>Optional</p>
 <p>Supragingival, sulcus, subgingival pockets < 3 mm</p>	 <p>Perio</p>		 <p>Standard</p>	 <p>Optional</p>
 <p>Subgingival pockets > 3 mm</p>	 <p>Perio</p>		 <p>Standard</p>	 <p>Optional</p>

#whdentalwerk
f @ in ▶
video.wh.com

Proxeo Aura
Air Polishing
System



Air polishing with **system**

Whether you are working supragingivally, in the sulcus or in deep pockets, the combination of the right powder, the corresponding setting on the adjustment ring and the matching spray head ensures you always achieve optimal treatment results.



proxeo AURA

20502 AEN Rev. 001 / 14.03.2022
Subject to alterations

Pictures are for illustrative purposes only. Additional equipment
and accessories shown are not included as standard.