

3M Science.
Applied to Life.™

Cement Selection Guide



Indications & Materials*

*Only primary cement recommendations
are shown. See respective instructions
for use for full list of indications.

Metal based
Zirconia
Glass Ceramics
Resin based

3M Deutschland GmbH, ESPE Platz, 82229 Saargölz, Germany



Veneers

- 3M™ RelyX™ Veneer
- 3M™ RelyX™ Universal
- 3M™ RelyX™ Universal
- 3M™ RelyX™ Veneer

Inlays/Onlays

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Veneer

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus

Bridges

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus

Crowns

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ RelyX™ Universal

Endodontic Posts

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ Ketac™ Cem Plus
- 3M™ Ketac™ Temp NE

Maryland Bridges

- RelyX™ Universal
- RelyX™ Unicem/Unicem 2

- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2
- 3M™ RelyX™ Universal
- 3M™ RelyX™ Temp NE
- 3M™ Ketac™ Cem Plus
- 3M™ RelyX™ Universal
- 3M™ RelyX™ Unicem/Unicem 2

On Implant Abutments

- 3M™ RelyX™ Universal

3M Science.
Applied to Life.™

Cement Selection Guide



Features & Benefits

3M, Ketac and RelyX are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
All other trademarks are owned by other companies. © 2021, 3M. All rights reserved. (02.2021)

3M™ RelyX™ Temp NE/E
The temporary cement.



- Reliable, temporary adhesion
- Quick, easy excess removal
- Easy removal of provisional and non-eugenol formulation

- Tack light-cure feature for "clean-up on demand"
- Easy to use automix delivery
- Sustained fluoride release
- Virtually no post-operative sensitivity



3M™ Ketac™ Cem Plus
Everyday easy.

- Saves time by eliminating etching, priming and bonding steps
- Virtually no post-operative sensitivity
- Strong and moisture-tolerant
- Endo tips make cementing posts easier and more reliable
- The world's most clinically proven self-adhesive resin cement



3M™ RelyX™ Unicem/3M™ RelyX™ Unicem 2
All around reliability.

- Customer-preferred shading system
- Extremely colour stable
- Matching try-in pastes
- Light-cure formula
- Easy handling and clean-up

3M™ RelyX™ Universal
Truly universal.

- Game-changing syringe
- Higher radiopacity
- Easy excess clean-up
- Improved bond strength
- 2-in-1: Adhesive and self-adhesive resin cement



3M™ RelyX™ Veneer
The veneer cement.